TML\MSH Microbiology Department	Policy #MI\ANTI\v23	Page 1 of 3
Policy & Procedure Manual		
Section: Antimicrobial Susceptibility Testing	Subject Title: Table of Contents	5
Manual		
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000	
Approved by: Laboratory Director	Revision Date: January 11, 2006	
	Review Date: April 21, 2005	

TABLE OF CONTENTS

Criteria for Susceptibility Testing
Verification of Antimicrobial Susceptibility Test Results and Confirmation of Organism
Identification
Reporting:
Urine - Gram Positive Susceptibility Reporting -1- Staphylococcus species, MRSA11
Urine - Gram Positive Susceptibility Reporting –2– Enterococcus, Streptococcus
Urine - Gram Negative Susceptibility Reporting
Enterics
Respiratory and Miscellaneous Non-Sterile Sites - Gram Positive Susceptibility Reporting -1-
Staphylococcus
Respiratory and Miscellaneous Non-Sterile Sites - Gram Positive Susceptibility Reporting – 2 –
Enterococcus, Streptococcus, Corynebacterim spp., Bacillus spp., viridans Streptococcus, Listeria spp.
Respiratory and Miscellaneous Non-Sterile Sites - Gram Negative Susceptibility Reporting -1-
Enterobacteriaceae, P. aeruginosa
Respiratory and Miscellaneous Non-Sterile Sites - Gram Negative Susceptibility Reporting – 2 –
Haemophilus species, Neisseria gonorrhoeae, S. maltophilia, B. cepacia, Pseudomonas species
(other than P. aeruginosa), fastidious gram-negative bacteria, non-fermenters, Neisseria
meningitidis
Spinal Fluids – Gram Positive Susceptibility Reporting
Spinal Fluids – Gram Negative Susceptibility Reporting –1– Enterobacteriaceae and Acinetobacter spp.,
Salmonella species
Spinal Fluids – Gram Negative Susceptibility Reporting –2– Pseudomonas aeruginosa, Pseudomonas spp.
(other than P. aeruginosa), fastidious gram-negative bacteria and non-fermenters, N. gonorrhoeae,
N. meningitidis Stenotrophomonas maltophilia, Burkholderia cepacia, Haemophilus species21
Blood and other Sterile Sites - Gram Positive Susceptibility Reporting -1- Staphylococcus species,
Enterococcus species

PROCEDURE MANUAL TORONTO MEDICAL LABORATORIES / MOUNT SINAI HOSPITAL MICROBIOLOGY DEPARTMENT

TML\MSH Microbiology Department	Policy #MI\ANTI\v23	Page 2 of 3
Policy & Procedure Manual		
Section: Antimicrobial Susceptibility Testing	Subject Title: Table of Contents	
Manual		

Blood and other Sterile Sites - Gram Positive Susceptibility Reporting – 2 - S. pneumoniae, viridans
Streptococcus, Streptococcus bovis,, S. milleri group, Group A, B, C, G Streptococcus, Listeria
species Corynebacterim species, Bacillus species
Blood and other Sterile Sites - Gram Negative Susceptibility Reporting -1 - Enterobacteriaceae and
Acinetobacter spp., Salmonella species including S. typhi
Blood and other Sterile Sites - Gram Negative Susceptibility Reporting -2 - Pseudomonas aeruginosa,
Pseudomonas species (other than P. aeruginosa), fastidious gram-negative bacteria, non-
fermenters, N. gonorrhoeae, N. meningitidis, Stenotrophomonas maltophilia, Burkholderia
cepacia, Haemophilus species
Antimicrobial related LIS Canned Messages
APPENDICES:
Disk Diffusion
Double Disk Diffusion for Erythromycin and Clindamycin on ß-haemolytic Streptococci and
Streptococcus pneumoniae
Double Disk Test for ESBL
beta-lactamase Testing
Oxacillin Screen Plate
Denka MRSA Screen
Serum Bacteriostatic & Bactericidal titres
Broth Macrodilution and Agar Dilution
Broth Microdilution MIC
E-Test
Vancomycin & High Level Aminoglycoside Testing
Vancomycin Screen Plate for <i>Staphylococcus</i>
Antimicrobial Abbreviations:
Disks
E-Test Strips
LIS

TML\MSH Microbiology Department Policy & Procedure Manual	Policy #MI\ANTI\v23	Page 3 of 3
Section: Antimicrobial Susceptibility Testing	Subject Title: Table of Contents	
Manual		

Oxacillin and Vancomycin Screen Recording Sheet for S. aureus	66
Quad Screen Recording Sheet for Enterococcus and CNST	67
VITEK	68
Record of Edited Revisions	72

TML\MSH Microbiology Department	Policy #MI\ANTI\01\v13 Page 1 c	
Policy & Procedure Manual		
Section: Antimicrobial Susceptibility Testing	Subject Title: Criteria for Susceptibility	
Manual	Testing	
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000	
Approved by: Laboratory Director	Revision Date: January 11, 2006)

CRITERIA FOR SUSCEPTIBILITY TESTING

I. <u>Introduction</u>

This section lists the susceptibility testing methods and required antimicrobials for each significant organism appropriate to the site of isolation.

II. <u>Reagents/Materials/Media</u>

Analytical Process - Bacteriology Reagents_Materials_Media List QPCMI10001

III. <u>Method</u>

- 1. Select significant organisms as per procedure manual of body sites.
- 2. Identify the selected isolate as per <u>Bacteria and Yeast Workup Manual</u>
- 3. For identical organisms, as defined in Bacteria and Yeast Workup Manual Minimal workup, isolated within 1 day from blood and sterile sites OR within 3 days from other sites do not require repeat susceptibility testing EXCEPT oxacillin and vancomycin screens for *Staphylococcus* and *Enterococcus*.
- 4. Refer susceptibility results back to like sites only and NEVER refer a sterile site to a non-sterile site.
- 5. Refer the susceptibility result to the previous cultures with the statement "Susceptibility testing not done. Please refer to No. (<u>HIS#</u> if present or <u>LIS#</u>) collected on <u>date</u> ".
- 6. Follow the table below as a guide for the appropriate method(s)/antimicrobial(s) to be set up.

Organisms	Site	Method	Antimicrobial(s)	
Aerobic Gram negative: Enterobacteriaceae, <i>Pseudomonas aeruginosa</i> and	All sites	Vitek ?	gns-623	
Acinetobacter species	If CPD=I or R, <i>E. coli</i> , <i>Klebsiella</i> species or <i>Proteus</i> species	KB ? add KB ?	FEP (for all PMH ONLY -non-CSF specimens - if CRO is I or R or <i>P.aeruginosa</i>) AMC, ASM, CRO, TAX, TAZ, CPD, FOX, ROX, TZP, FEP	
<i>E. coli</i> O157	Enteric sites	Not tested		
All Salmonella species	Blood and Sterile sites	add KB ?	NA	
PROCEDURE MANUAL				

TML\MSH Microbiology Department	Policy # MI\ANTI\01\v13	Page 2 of 4
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual		

Organisms Salmonella species other than S. typhi	Site Enterics from TG and TW	Method Vitek ?	Antimicrobial(s) gns-623
	Enterics from other hospitals	Not tested	
Stenotrophomonas maltophilia	All sites	KB ?	SXT, LVX
Burkholderia cepacia	All sites	KB ?	SXT, TAZ, MEM
Pseudomonas aeruginosa, Acinetobacter species Stenotrophomonas maltophilia or Burkholderia cepacia that is resistant to all routinely tested antimicrobials	All sites	e-test ?	PO
<i>Pseudomonas</i> sp., <i>Aeromonas</i> sp., <i>Plesiomonas</i> sp. and other afermenters	All sites	Not tested	
Mucoid <i>P. aeruginosa,</i> <i>Enterobacteriaceae</i> not growing on Vitek	All sites	KB ?	AMP, KZ, CRO, CIP, SXT, GM, TZP, TOB, CAZ, IMI, CPD
Haemophilus species	All sites	beta-lactamase	
	Blood and Sterile sites	beta-lactamase +	CRO CID AMD
Moraxella catarrhalis	All sites	Not tested	CKO, CIF, AMIF
Neisseria gonorrhoeae	All sites	Not tested	
Neisseria meningitides	All sites	Not tested	
Other fastidious Gram negatives (e.g. <i>Pasteurella</i> species)	All sites	Not tested	
Campylobacter species	All sites	Not tested	
Staphylococcus aureus	All sites	Vitek ? + Screen plate ?	gps-105 OX, VA
MRSA – add e-test MU until further notice. (KB-MUP pending validation)	MRSA	a ? add e-test ?	MU
	If MRSA from MRSA Screen Test and TE and SXT=R	add e-test ?	FA

TML\MSH Microbiology Department	Policy # MI\ANTI\01\v13	Page 3 of 4
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual		

Organisms	Site	Method	Antimicrobial(s)
Coagulase-negative	Urine	Not tested	
Staphylococcus	All other sites (test on	Vitek ?	gps-105
	pooled isolates if more	+	5p5 100
	than one morphotype of	Screen plate ?	VA
	CNST)		
Micrococcus species	All sites	Not tested	
Enterococcus species	Urines	Vitek ?	gps-105
		Screen plate ?	VA
	Blood & Sterile sites	<u>KB</u> ?	AMP
		+	
		Screen plate ?	High level GM and
		+ beta-lactamase	STREP, VA
	All other sites	KB ?	AMP
		+	
		<u>Screen plate</u> ?	<u>VA</u>
	If VA=R, <i>E. faecalis</i> or <i>E. faecium</i>	add VRE MIC panel	
Streptococcus	Blood & Sterile sites	e-test ?	P, CRO
pneumoniae		+ Double disc KP 9	CC F
			CC, E
		KB ?	LVX, VA
	All other sites	Double disc KB ?	CC, E
		н КВ ?	OX, LVX, VA
		If OX= \mathbf{R} , then e-test ?	P, CRO
Group A, B, C, G	? Blood and Sterile sites	Double disc KB ?	CC, E
Streptococcus		+ VD 9	D VA
	? Urine, on request	<u> </u>	$\frac{\mathbf{r}, \mathbf{v}\mathbf{A}}{\mathbf{L}\mathbf{V}\mathbf{X}}$
	ONLY		
	? Vaginal GBS screens,	Double disc KB ?	CC, E
	on request ONLY	+	
		KB ?	VA
	? Other sites, on request	Double disc KB ?	CC, E
	UNL Y	⁺ KB ?	LVX, VA

TML\MSH Microbiology Department	Policy # MI\ANTI\01\v13	Page 4 of 4
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual		

Organisms	Site	Method	Antimicrobial(s)
Streptococcus bovis	Blood & Sterile sites	e-test ?	CRO, P, VA
	Mixed cultures	Not tested	
	Other sites	Not tested	
viridans Streptococci	Blood & Sterile sites	e-test for pure cultures?	CRO, P, VA
	Mixed cultures	Not tested	
	Other sites	Not tested	
Streptococcus milleri group	Blood & Sterile sites	e-test ?	CRO, P, VA
<u> </u>	Urine	KB ?	P, LVX
	Other sites	KB ?	CC, E, LVX, P, VA
Listeria species	All sites	Not tested	
Corynebacterium species	All sites	Not tested	
Bacillus species	All sites	Not tested	
Nocardia species	All sites	Not tested	(Sent to PHL on special
Anaerobes	All sites	Not tested	request) (Sent to PHL on special request)

TML\MSH Microbiology Department	Policy #MI\ANTI\02\01\v01	Page 1 of 3	
Section: Antimicrobial Susceptibility Testing	Subject Title: Verification of Antimicrobia		
Manual	Susceptibility Test Results and		
	Confirmation of Organism Ida	ntification	
	Commination of Organism fue	initication	
Issued by: LABORATORY MANAGER	Original Date: November 21, 20	05	
Issued by: LABORATORY MANAGER Approved by: Laboratory Director	Original Date: November 21, 20 Revision Date:	05	

I. <u>Introduction</u>

This section describes the occasions where the drugs tested against isolates showed phenotype that:

- 1. have never been documented
- 2. are uncommon, and/or
- 3. represent results that could easily occur from technical errors and which may have significant clinical consequences.

II. <u>Reagents/Materials/Media</u>

Analytical Process - Bacteriology Reagents_Materials_Media List QPCMI10001

III. <u>Procedure</u>

When any of the listed results in the <u>TABLE 1</u> below occurs, verify the result as follows:

- 1. Check purity plate.
- 2. Check previous reports on the patient.
- 3. Confirm the identification of the isolate from the original isolation medium.
- 4. Repeat susceptibility test to confirm result. Use an alternative method if applicable.
- For isolates that show results other than susceptible for those antimicrobial agents for which only susceptible interpretive criteria are provided by NCCLS guidelines M100-S15 (listed as "not S" above) an for staphylococci with vancomycin – I or R results:
 - i Confirm the organism identification
 - ii. Confirm the antimicrobial susceptibility test results
 - iii. Freeze the isolate
 - iv. Send the isolate to PHL for confirmation.
- 6. If the result is confirmed, notify the Charge Technologist.
- 7. The Charge Technologist confirms the result and notifies the Microbiologist.
- 8. The Microbiologist further confirms the result and notifies the Infection Control Practitioner.

For results marked with *, LIS reflex rules will automatically report these as R; repeat susceptibility testing is not required if the purity and organism identification is confirmed.

TML\MSH Microbiology Department Policy & Procedure Manual Policy # MI\ANTI\02\01\v01 Page

Page 2 of 3

Antimicrobial Susceptibility Testing Manual

Organism or Group	Uncommon results
Any organism	Resistant to all agents routinely tested
Gram-negative organisms	
Enterobacteriaceae	Imipenem, Meropenem – I or R
Citrobacter freundii	Ampicillin, Cefazolin – S*
Enterobacter species	
Serratia marcescens	
Klebsiella species	Ampicillin – S*
Proteous vulgaris	
Providencia species	
Escherichia coli, Klebsiella species,	Cefpodoxime – Vitek=I or R; KB=S
Proteus species	
Escherichia coli, Klebsiella species,	KB-ESBL panel with reduction in zone
Proteus species	of inhibition instead of potentiation or
	no change.
Stenotrophomonas maltophilia	Imipenem, Meropenem – S
Haemophilus influenzae	Aztreonam – not S
	Imipenem, Meropenem – not S
	3 rd generation cephalosporin – not S
	Fluoroqinolone – not S
Neisseria gonorrhoeae	3 rd generation cephalosporin – R
Gram-positive organisms	
Enterococcus faecalis	Ampicillin or Penicillin – R
	Daptomycin – not S
	Quinupristin-Dalfopristin – S
	Linezolid – I or R
Enterococcus faecium	Daptomycin – not S
	Linezolid – I or R
Staphylococcus aureus	Daptomycin – not S
	Linezolid – not S
	Quinupristin-Dalfopristin – I or R
	Vancomycin – I or R
Coagulase-negative Staphylococcus	Daptomycin – not S
	Linezolid – not S
	Vancomycin – I or R
Streptococcus pneumoniae	3^{iu} generation cephalosporin – R
	Fluoroquinolone - R
	Linezolid – not S
	Vancomycin – not S

TML\MSH Microbiology Department	Policy # MI\ANTI\02\01\v01	Page 3 of 3
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual		

Organism or Group	Uncommon results
Gram-positive organisms (cont'd)	
beta-haemolytic Streptococcus	Ampicillin or Penicillin – not S
	3 rd generation cephalosporin – not S
	Daptomycin – not S
	Linezolid – not S
	Vancomycin – not S
viridans Streptococcus	Daptomycin – not S
	Linezolid – not S
	Vancomycin – not S

IV. <u>Reference</u>

Suggestions for Verification of Antimicrobial susceptibility Test Results and Confirmation of Organism identification in Table 8 of Clinical and Laboratory Standards Institute (CLSI) Document - Performance Standards for Antimicrobial Susceptibility Testing M100-S15, Vol. 25 No.1, 2005.

TML\MSH Microbiology Department	Policy #MI\ANTI\03\01\v10	Page 1 of 1			
Policy & Procedure Manual					
Section: Antimicrobial Susceptibility Testing	Subject Title: Urine - Gram Positive				
Manual	Susceptibility Reporting – 1				
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000				
Approved by: Laboratory Director	Revision Date: November 21, 2005				
	Review Date: April 21, 2005				
Site: Uning Crom Disting Suscentibility Departing 1 Stanbulgeorgy apoping MDSA					

Site: Urine – Gram Positive Susceptibility Reporting – 1 – Staphylococcus species, MRSA

	Staphylococcus aureus				MRSA	
Antimicrobial Agent	MSH ^a	UHN ^b	CHC ^c	MSH ^a	UHN ^b	CHC
Ampicillin	X ¹	X^1	X ¹	X ^{1,5}	X ^{1, 5}	X ^{1, 5}
Cefazolin	X ²	X ²	X ²	X ^{2, 5}	X ^{2, 5}	X ^{2, 5}
Cloxacillin	X ²	X ²	X ²	X ^{2, 5}	X ^{2, 5}	X ^{2, 5}
Doxycycline				X ^{3, 6}	X ^{3, 6}	X ^{3, 6}
Fusidic Acid				X ^{6, 7}	X ^{6, 7}	X ^{6, 7}
Mupirocin				X ⁶	X ⁶	X ⁶
Nitrofurantoin	Х	Х	Х	X ⁵	X ⁵	X ⁵
Rifampin				X ⁶	X ⁶	X^6
Trimethoprim/Sulfa	X	Х	Х	Х	Х	Х
Vancomycin	X ⁴	X ⁴	X ⁴	X ⁵	X ⁵	X ⁵

¹ Base on Penicillin or beta-lactamase result

² Base on Oxacillin result

³ Adults only (>13 y); base on Tetracycline result

⁴ Report if patient is allergic to Penicillin OR if *Staphylococcus* species is resistant to **All** other antimicrobial agents.

⁵ DO NOT report if isolated from Infection Control Screening test

⁶ Report ONLY if isolated from Infection Control Screening test. Include Isolate Comment "Susceptibility results are provided to guide decolonization therapy; if decolonization therapy is being considered, please consult your infection control team."

⁷ Report only if resistant to Mupirocin.

^a includes MSH, PMH, Baycrest, TRI, CAMH

^b includes TG, TW, Bridgepoint, Grace

^c includes CHC, Ajax/Pickering

If all antimicrobial agents are resistant, inform the Microbiologist on-call.

Note: *S. saprophyticus* and **coagulase-negative***-Staphylococcus* - **DO NOT** report susceptibilities. Report with Isolate Comment - "Susceptibility testing of this organism is not warranted because infections respond to concentrations achieved in urine of antimicrobial agents commonly used to treat acute, uncomplicated urinary tract infections e.g. nitrofurantoin, trimethoprim/sulfa or fluoroquinolones.

TML\MSH Micr	ML\MSH Microbiology Department		Policy #MI\ANTI\03\02\v10	Page 1 of 1]		
Policy & Procedure Manual							
Section: Antimic	crobial Su	sceptibility	y Testing	Subject Title: Urine - Gram Posi	Subject Title: Urine - Gram Positive		
Manual				Susceptibility Repo	orting – 2		
Issued by: LAB	ORATOR	Y MANA	GER	Original Date: January 10, 2000			
Approved by: La	aboratory D	Director		Revision Date: November 21, 200)5		
				Review Date: April 21, 2005			
Urine – Gram F	Positive Su	sceptibilit	y Reporting	g – 2 – Enterococcus, Streptococcu	5		
	Enter	<i>ococcus</i> sp	ecies ¹¹	Group A, B, C, G Streptococcus	Streptoco mille	occus ri	
Antimicrobial Agent	MOLIA	TTTT	aua	All sites			
Agent	MSH"	UHN	CHC	Routinely not tested. See below ⁸ .	All sit	es	
				For special request:			
Ampicillin	Х	Х	Х				
Ciprofloxacin	$X^{4,6}$	$X^{5,6}$	X^6				
Clindamycin							
Erythromycin							
Levofloxacin				X^{10}	X ⁶		
Linezolid	X^2	X^2	X^2				
Nitrofurantoin	Х	Х	Х				
Penicillin G					X		
Synercid	X ²	X^2	X^2				
Tetracycline	X ¹	X ¹	X ¹				
Vancomycin	X ^{3, 9}	X ^{3, 9}	X ^{3, 9}				

¹ Adults only (>13 y) ² If Vancomycin, Ciprofloxacin, Tetracycline, Nitrofurantoin and Ampicillin are R except for *E. gallinarum* and *E. casseliflavus*.

³ Test but **DO NOT** report unless Vancomycin R or Enterococcus resistant to **All** other antimicrobial agents ⁴ TRI only

⁵ Test and report on request only

⁶ Adults only (>18y) ⁷ Always report as sensitive

⁸ Report "This organism is intrinsically susceptible to penicillin. If treatment is required and this patient cannot be treated with penicillin, please contact the Microbiology Department within 48 hours to request sensitivity testing."

⁹ E. gallinarum and E. casseliflavus, report as \mathbf{R} with the statement "This organism always has intrinsic nontransmissible resistance to vancomycin. The patient does not require isolation." ¹⁰Report with additional isolate comment "Susceptibility completed as requested" (do not remove original

comments). If Levofloxacin is R or patient is <18y, consult the Microbiologist.

¹¹ DO NOT report susceptibilities if isolated from Infection Control Screening test.

^a includes MSH, PMH, Baycrest, TRI, CAMH

^b includes TG, TW, Bridgepoint, Grace

^c includes CHC, Ajax/Pickering

Note: If all antimicrobial agents are resistant, inform the Microbiologist on-call.

TML\MSH Microbiology Department Policy & Procedure Manual	Policy #MI\ANTI\03\03\v10	Page 1 of 1		
Section: Antimicrobial Susceptibility Testing	Subject Title: Urine - Gram Negative			
Manual	Susceptibility Reporting			
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000			
Approved by: Laboratory Director	Revision Date: November 21, 2005			
	Review Date: April 21, 2005			

Antimicrobial Agent	Entero	terobacteriaceae and Pseudomonas		S. maltophilia	Burkholderia			
	Acinetobacter spp.			aeruginosa				cepacia
	MSH ^a	UHN ^b	CHC ^c	MSH ^a	UHN ^b	CHC ^c	All Sites	All Sites
Amikacin	X^{14}	X ¹⁴	X^{14}	X ¹⁴	X^{14}	X^{14}		
Ampicillin	X', 13	X', 13	X', 13	\mathbf{R}^2				
Cefazolin	$X^{8,10,13}$	$X^{8, 10, 13}$	X ^{8, 10, 13}	\mathbf{R}^2				
Cefepime	$X^{11,13}$			X^{12}				
Ceftazidime				Х	Х	Х		Х
Ceftriaxone	X ^{6, 9, 13}	X ^{6, 9, 13}	X ^{6, 9, 13}					
Ciprofloxacin	X ¹	$X^{1,5}$		Xi	$X^{1,5}$			
Colistin	X ¹⁵	X ¹⁵	X ¹⁵	X ¹⁵	X^{15}	X ¹⁵	X ¹⁵	X ¹⁵
Gentamicin	Х	Х	Х	Х	Х	Х		
Imipenem	X^4	X^4	X^4	X^4	X^4	X^4		
Levofloxacin							X	
Meropenem								Х
Nitrofurantoin	X ¹⁰	X ¹⁰	X ¹⁰	R^2				
Ofloxacin		X ^{1,5}	X ^{1,5}		$X^{1,5}$	X ^{1,5}		
Piperacillin/Tazobactam	$X^{9, 10}$	X9		Х	Х	Х		
Trimethoprim/Sulfa	Х	Х	Х	\mathbf{R}^2			Х	Х
Tobramycin	Х	X ³	X ³	Х	X^3	X ³		
¹ Adults only (>18 y)		2 R = A	Always repo	ort as resi	stant		³ Report if C	Genta R

¹ Adults only (>18 y) 2 R = Always report as resistant ⁴ Report if R **OR** if R to **All** other Antimicrobial Agents **OR** if only aminoglycoside is S

⁵ Report Ofloxacin (base on Ciprofloxacin result) for Enterobacteriaceae on all Bridgepoint, CHC and Ajax/Pickering patients. **DO NOT** report Ciprofloxacin.

⁶ Report ceftriaxone only for *Acinetobacter* spp. or if *Enterobacteriaceae* is intermediate or resistant to ceftriaxone

⁷ Klebsiella spp., Enterobacter spp., Acinetobacter spp., H. alvei, Citrobacter spp., Pantoea agglonerans, Proteus vulgaris, Proteus penneri, & Serratia species always report Amp as R.

⁸ Enterobacter spp., Citrobacter spp., Pantoea agglonerans, Proteus vulgaris, Proteus penneri, Acinetobacter spp. & Serratia species always report Cefazolin and Cephalothin as R.

⁹ Cedecea spp., Citrobacter spp., Enterobacter spp., Hafnia spp., Morganella morganii, Pantoea agglonerans, Proteus penneri, Proteus vulgaris, Providencia species, Seratia species, if S, report with comment "Resistance to extended-spectrum penicillins, beta-lactam/beta-lactamase inhibitor combinations, and cephalosporins may develop during therapy with these agents. For serious infections, these agents should be avoided and consultation with a medical microbiologist or infectious disease physician is strongly recommended."

¹⁰ Do not report for Salmonella species.

¹¹ Report on PMH patients only if Ceftriaxone is I or R.

¹² PMH only

¹³ For *Enterobacteriaceae* if any one of cefotaxime/ceftriaxone or ceftazidime=R, report all as R

¹⁴Report if both Gentamicin and Tobramycin are R.

¹⁵For isolates other than Enterobacteriaceae, report if R to all other drugs; base on Polymyxin B.

^a includes MSH, PMH, Baycrest, TRI, CAMH ^b includes TG, TW, Bridgepoint, Grace ^c includes CHC, Ajax/Pickering

Note: Pseudomonas species (other than P. aeruginosa), fastidious gram-negative bacteria & non-fermenters - DO

NOT report susceptibility result. Report with ISOLATE comment "In vitro susceptibility testing for this organism is not routinely performed and/or is unreliable. If advice on antimicrobial therapy is required, please contact the Medical Microbiologist".

If all antimicrobial agents are resistant, inform the Microbiologist on-call.

PROCEDURE MANUAL

TML\MSH Microbiology Department	Policy #MI\ANTI\03\04\v10	Page 1 of 1	
Policy & Procedure Manual			
Section: Antimicrobial Susceptibility Testing	Subject Title: Enterics Susceptibility		
Manual	Reporting		
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000		
Approved by: Laboratory Director	Revision Date: November 21, 20	005	
	Review Date: April 21, 2005		

Enterics

Antimicrobial	Shig	gella spe	ecies	Salm	o <i>nella</i> sp	oecies	Saln	ionella i	typhi	Vibr	io chole	erae ³
Agent	MSH ^a	UHN ^b	CHC ^c	MSH ^a	UHN ^b	CHC ^c	MSH ^a	UHN ^b	CHC ^c	MSH ^a	UHN⁵	CHC ^c
Ampicillin	X	Х	Х		X^2		Х	Х	Х	Х	Х	Х
Cefepime							X°					
Ceftriaxone							Х	Х	Х			
Ciprofloxacin	X	X	X		$X^{1,2}$		X	X	X	X	X	X
Trimethoprim/Sulfa	X	Х	Х		X^2		Х	Х	X	Х	Х	Х
Tetracycline										X^4	X^4	X^4

¹ Adults only (>18 y)
² For TG and TW patients, test but **DO NOT** report.
³ For all Vibrio species other than *V. cholerae*, **DO NOT** test or report.

⁴ Adults only (>13 y)

⁵ Report on PMH patients only if Ceftriaxone is I or R, base on Ceftriaxone result.

^a includes MSH, PMH, Baycrest, TRI, CAMH

^b includes TG, TW, Bridgepoint, Grace

^c includes CHC, Ajax/Pickering

Note: E. coli O157, Campylobacter spp., and Yersinia spp. - DO NOT report susceptibility result. Report with ISOLATE comment "In vitro susceptibility testing for this organism is not routinely performed and/or is unreliable. If advice on antimicrobial therapy is required, please contact the Medical Microbiologist".

If all antimicrobial agents are resistant, inform the Microbiologist on-call.

TML\MSH Microbiology Department	Policy #MI\ANTI\03\05\v10	Page 1 of 1			
Policy & Procedure Manual					
Section: Antimicrobial Susceptibility Testing	Subject Title: Respiratory and	Miscellaneous			
Manual	Non-Sterile Sites - Gram Positiv				
	Susceptibility Reporting – 1				
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000				
Approved by: Laboratory Director	Revision Date: November 21, 20	005			
	Review Date: April 21, 2005				

Respiratory and Miscellaneous Non-Sterile Sites - Gram Positive Susceptibility Reporting -1 -**Staphylococcus**

Antimicrobial Agent	Staph	ylococcus sp	oecies		MRSA	
Antimier obiar Algent	MSH ^a	UHN ^b	CHC ^c	MSH ^a	UHN⁵	CHC ^c
Cefazolin	X^2	X^2	X ²	X ^{2, 6}	X ^{2, 6}	X ^{2, 6}
Clindamycin	X^4	X^4	X ⁴	X ^{4, 6}	X ^{4, 6}	X ^{4, 6}
Cloxacillin	X^2	X^2	X ²	X ^{2, 6}	X ^{2, 6}	X ^{2, 6}
Doxycycline				X ^{3, 5}	X ^{3, 5}	X ^{3, 5}
Erythromycin	X^4	X^4	X ⁴	X ^{4, 6}	X ^{4, 6}	X ^{4, 6}
Fusidic Acid				X ^{3, 7}	X ^{3, 7}	X ^{3, 7}
Mupirocin				X ³	X ³	X ³
Rifampin				X ³	X^3	X ³
TMP/SMX	Х	Х	X	Х	Х	Х
Vancomycin	X^1	X ¹	X ¹	X ⁶	X^6	X ⁶

¹ Report if Oxacillin R

² Base on Oxacillin result

³ Report ONLY if isolated from Infection Control Screen. Include Isolate Comment "Susceptibility results are provided to guide decolonization therapy; if decolonization therapy is being considered, please consult your infection control team."

 4 Clinda R = Erythro R

⁵ Adults only (>13yrs); base on Tetracycline result

⁶ DO NOT report if isolated from Infection Control Screen.

⁷ Report only if resistant to Mupirocin.

^a includes MSH, PMH, Baycrest, TRI, CAMH

^b includes TG, TW, Bridgepoint, Grace ^c includes CHC, Ajax/Pickering

Note: For organisms isolated from ears and eyes sources and susceptibility result is reported, add comment "These susceptibility testing results are based on guidelines for systemic antimicrobial agents and may not accurately represent activity of topical agents."

If all antimicrobial agents are resistant, inform the Microbiologist on-call.

TML\MSH Microbiology Department	Policy #MI\ANTI\03\06\v10	Page 1 of 1				
Policy & Procedure Manual						
Section: Antimicrobial Susceptibility Testing Manual	Subject Title: Respiratory and Miscellaneous Non-Sterile Sites -					
	Gram Positive Susceptibility Reporting – 2					
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000					
Approved by: Laboratory Director	Revision Date: November 21, 2005					
	Review Date: April 21, 2005					

Respiratory and Miscellaneous Non-Sterile Sites - Gram Positive Susceptibility Reporting-2-Enterococcus, Streptococcus, Corynebacterim spp., Bacillus spp., viridans Streptococcus, Listeria spp.

Antimicrobial Agent	En	Enterococcus ¹		S. pneumoniae		iae	Group A, B, C, G Streptococcus	S. milleri
	MSH ^a	UHN ^b	CHC ^c	MSH ^a	UHN [♭]	CHC ^c	All Sites	All Sites
Ampicillin	X	Х	X				Routinely not tested. See	
Cefepime				X ¹⁴			below ⁷ .	
Ceftriaxone				X°	X	X°	For special request:	
Clindamycin				X^4	X^4	X^4	X ^{4, 12}	Х
Erythromycin				X ^{4, 5}	X ^{4, 5}	X ^{4, 5}	X ^{4, 12}	Х
Levofloxacin				X ^{10, 11}		X ⁱⁱ	$X^{2, 11, 12, 15}$	Х
Linezolid	X ¹³	X ¹³	X ¹³					
Moxifloxacin				X ^{9, 11}	X ^{9, 11}	_		
Penicillin G				X ⁸	X ⁸	X ⁸		Х
Synercid	X ¹³	X ¹³	X ¹³					
Vancomycin	X ³	X ³	X ³	X°	X ⁶	X	$X^{12,16}$	

¹ DO NOT report if isolated from Infection Control Screen.

² DO NOT report on GBS Screen or vaginal swab

³ E. gallinarum, and E. casseliflavus, report as \mathbf{R} with the statement "This organism always has intrinsic non-transmissible resistance to vancomycin. The patient does not require isolation."

⁴ Report as R if D-zone is present. Report clindamycin as R if erythromycin is R

⁵ Report Erythromycin for respiratory specimens only

⁶ Report only if Pen I or R

⁷ Report "This organism is intrinsically susceptible to penicillin. If treatment is required and this patient cannot be treated with penicillin, please contact the Microbiology Department within 48 hours to request sensitivity testing."

⁸ Base on Oxacillin result if S; base on Penicillin Etest if Oxacillin is R

⁹ Base on Levofloxacin result. Report on MSH, PMH, TG and TW patients only.

¹⁰ DO NOT report on MSH and PMH patients.

¹¹Adults only (>18 yrs)

¹²Report with additional isolate comment "Susceptibility completed as requested" (do not remove original comments).

¹³ If Vancomycin and Ampicillin are R except for *E. gallinarum* and *E. casseliflavus*.
 ¹⁴ Report on PMH patients only if Cetriaxone is I or R, based on Ceftriaxone result.
 ¹⁵ If Levofloxacin is R or patient is <18y, consult the Microbiologist.

¹⁶Report only if either Clindamycin or Erythromycin are I or R. ^a includes MSH, PMH, Baycrest, TRI, CAMH

^b includes TG, TW, Bridgepoint, Grace

^b includes TG, TW, Bridgepoint, Grace ^c includes CHC, Ajax/Pickering Note: *Listeria* species – DO NOT report susceptibility result. Report with ISOLATE comment –"Routine in vitro

susceptibility testing of this organism is unreliable. *Listeria* spp. should be considered resistant to all cephalosporins. The recommended regimen for therapy is ampicillin. If additional advice on antimicrobial therapy is required, please contact the Medical Microbiologist." Corynebacterim species, Bacillus species, viridans Streptococcus - DO NOT report susceptibility result. Report with

ISOLATE comment "In vitro susceptibility testing for this organism is not routinely performed and/or is unreliable. If advice on antimicrobial therapy is required, please contact the Medical Microbiologist".

For organisms isolated from ears and eyes sources and susceptibility result is reported, add comment "These susceptibility testing results are based on guidelines for systemic antimicrobial agents and may not accurately represent activity of topical agents."

If all antimicrobial agents are resistant, inform the Microbiologist on-call.

TML\MSH Microbiology Department	Policy #MI\ANTI\03\07\v10	Page 1 of 1	
Policy & Procedure Manual			
Section: Antimicrobial Susceptibility Testing Manual	Subject Title: Respiratory and Miscel	laneous Non-sterile	
	Sites - Gram Negative Susceptibility		
	Reporting – 1		
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000		
Approved by: Laboratory Director	Revision Date: November 21, 2005		
	Review Date: April 21, 2005		

Respiratory and Miscellaneous Non-Sterile Sites - Gram Negative Susceptibility Reporting – 1 – Enterobacteriaceae, *P. aeruginosa*

Antimicrobial Agent	Enter Aci	robacteriacea inetobacter s	e and pp.	Pseudo	omonas aeru	ginosa
	MSH ^a	UHN⁰	CHC ^c	MSH ^a	UHN⁰	CHC ^c
Amikacin	X^{14}	X ¹⁴	X^{14}	X^{14}	X^{i_4}	X^{14}
Ampicillin	$X^{2, 12, 13}$	$X^{2, 12, 13}$	$X^{2, 12, 13}$	R		
Cefazolin	$X^{3,12,13}$	X ^{3,12,13}	X ^{3,12,13}	R		
Cefepime	$X^{9;12}$			X^{10}		
Ceftazidime	$X^{4;12}$	$X^{4, 12}$	$X^{4,12}$	Х	Х	Х
Ceftriaxone	X ^{7,12}	X ^{7, 12}	X ^{7,12}			
Cefuroxime		$X^{8, 12}$	$X^{8, 12}$			
Ciprofloxacin	X ⁱ	X	X	X	X ⁱ	X ⁱ
Colistin	X ¹⁵	X ¹⁵	X ¹⁵	X ¹⁵	X ¹⁵	X ¹⁵
Gentamicin	Х	X	Х	Х	Х	Х
Imipenem	X ³	X ⁵	X ⁵	X ⁵	X	X ⁵
Piperacillin/Tazobactam	X ^{7, 11}	X^7		Х	Х	Х
Trimethoprim/Sulfa	Х	Х	Х	R		
Tobramycin	Х	X ⁶	X	Х	Х	Х

¹ Adults only (>18 y)

² Always report *Klebsiella* spp., *Enterobacter* spp. *Acinetobacter* spp., *H. alvei, Proteus vulgaris, Proteus penneri, Citrobacter* spp. and *Serratia* species as R

³ Always report Enterobacter spp., Citrobacter spp., Proteus vulgaris, Proteus penneri, Acinetobacter spp. and Serratia species as R.

⁴ Report only if R. For *Enterobacteriaceae* if cefotaxime/ceftriaxone or ceftazidime R, report both as R

⁵ Report if **R** OR if R to All other Antimicrobial Agents OR if only aminoglycoside is S

⁶ Report if Genta is R

⁷ *Cedecea* spp., *Citrobacter* spp., *Enterobacter* spp., *Hafnia* spp., *Morganella morganii, Pantoea agglonerans, Proteus penneri, Proteus vulgaris, Providencia* species, *Seratia* species, if S, report with comment "Resistance to extended-spectrum penicillins, beta-lactam/beta-lactamase inhibitor combinations, and cephalosporins may develop during therapy with these agents. For serious infections, these agents should be avoided and consultation with a medical microbiologist or infectious disease physician is strongly recommended."

⁸ Report on *Enterobacteriaceae*; Do not report for *Acinetobacter* spp. (Base on KB result) on respiratory specimens only.

⁹ Report on PMH patients only if Ceftriaxone is I or R.

¹⁰PMH only

¹¹Do not report for *Salmonella* species.

¹²For *Enterobacteriaceae* if any one of cefotaxime/ceftriaxone or ceftazidime=R, report all as R

¹³Report as I or R if Cefuroxime is I or R.

¹⁴Report if both Gentamicin and Tobramycin are R.

¹⁵For Acinetobacter species and *Pseudomonas aeruginosa*, report if R to all other drugs; base on Polymyxin B.

^a includes MSH, PMH, Baycrest, TRI, CAMH ^b includes TG, TW, Bridgepoint, Grace ^c includes CHC, Ajax/Pickering Note: *Pseudomonas* species (other than *P. aeruginosa*), fastidious gram-negative bacteria & non –fermenters - DO NOT

report susceptibility result. Report with ISOLATE comment "In vitro susceptibility testing for this organism is not routinely performed and/or is unreliable. If advice on antimicrobial therapy is required, please contact the Medical Microbiologist".

For organisms isolated from **ears and eyes sources** and susceptibility result is reported, add comment "These susceptibility testing results are based on guidelines for systemic antimicrobial agents and may not accurately represent activity of topical agents."

If all antimicrobial agents are resistant, inform the Microbiologist on-call.

PROCEDURE MANUAL

TML\MSH Microbiology Department Policy & Procedure Manual	Policy #MI\ANTI\03\08\v10	Page 1 of 1			
Section: Antimicrobial Susceptibility Testing	Subject Title: Respiratory and	Miscellaneous			
Manual	Non-Sterile Sites - Gram Negative				
	Susceptibility Rep	porting – 2			
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000				
Approved by: Laboratory Director	Revision Date: November 21, 20)05			
	Review Date: April 21, 2005				

Respiratory and Miscellaneous Non-Sterile Sites - Gram Negative Susceptibility Reporting – 2 – *Haemophilus* species, *M. catarrhalis, Neisseria gonorrhoeae, S. maltophilia, B. cepacia, Pseudomonas* species (other than *P. aeruginosa*), fastidious gram-negative bacteria, non-fermenters, *Neisseria meningitidis*

Antimicrobial Agent	Haemophilus species	S. maltophilia	Burkholderia cepacia
	All sites	All sites	All sites
Beta-lactamase	X^2		
Ceftazidime			Х
Levofloxacin		X ¹	
Meropenem			Х
Trimethoprim/Sulfa		Х	Х

¹ Adults only (>18 y)

² If beta-lactamase is negative, add comment "beta-lactamase negative result suggests susceptible to ampicillin."

If beta-lactamase is positive, add comment "beta-lactamase positive result suggests resistance to ampicillin but generally susceptible to amoxicillin-clavulanic acid and cefuroxime."

^a includes MSH, PMH, Baycrest, TRI, CAMH

^b includes TG, TW, Bridgepoint, Grace

^c includes CHC, Ajax/Pickering

Note: *Pseudomonas* species (other than *P. aeruginosa*), fastidious gram-negative bacteria, nonfermenters, *N. gonorrohoeae*, *N. meningitidis* - DO NOT report susceptibility result. Report with ISOLATE comment: "In vitro susceptibility testing for this organism is not routinely performed and/or is unreliable. If advice on antimicrobial therapy is required, please contact the Medical Microbiologist".

For *M. catarrhalis* - DO NOT report susceptibility result. Report with ISOLATE comment: "The majority of *Moraxella catarrhalis* are resistant to ampicillin. In vitro susceptibility testing for this organism is not routinely performed and/or is unreliable. If advice on antimicrobial therapy is required, please contact the Medical Microbiologist".

For organisms isolated from **ears and eyes sources** and susceptibility result is reported, add comment "These susceptibility testing results are based on guidelines for systemic antimicrobial agents and may not accurately represent activity of topical agents."

If all antimicrobial agents are resistant, inform the Microbiologist on-call.

TML\MSH Microbiology Department Policy & Procedure Manual	Policy #MI\ANTI\03\09\v10	Page 1 of 1			
Section: Antimicrobial Susceptibility Testing	Subject Title: Spinal Fluids – Gran	m Positive			
Manual	Susceptibility Reporting				
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000				
Approved by: Laboratory Director	Revision Date: November 21, 2005				
	Review Date: April 21, 2005				

Spinal Fluids – Gram Positive Susceptibility Reporting

Antimicrobial	Stapi	hyloco	ccus	En	terococi	cus	Strep.	viridans <i>Strep</i> .	S. milleri	Group A,B,C,G
A4	S	species			species		pneumoniae	S. Dovis	group	Sheptococcus
Agent	MSH ^a	UH№	CHC	MSH ^a	UHN⁵	CHC ^c	All Sites	All Sites	All Sites	All Sites
Ampicillin				Х	Х	Х				
Ceftriaxone							X ^{8, 4}	X ^{8, 4}	$X^{8, 4}$	
Cloxacillin	Х	Х	Х							
HLGR ³				Х	Х	Х				
HLSR ³				X^2	X^2	X^2				
Linezolid				X'	Χ′	Χ'				
Penicillin							X ⁸	X ⁸	X ⁸	Х
Piperacillin/Tazo				X ⁵						
Synercid				X'	Χ'	Χ'				
Trimethoprim/Sulfa	Х	Х	Х							
Tobramycin										_
Vancomycin	X	X	Х	X ⁶	X ⁶	X ⁶	X^4	X ⁴	$X^{8, 4}$	X^4

¹ Report if Oxacillin R

² Report only if requested.

³ HLGR = High Level Gentamicin Resistant; HLSR = High Level Streptomycin Resistant. Report base on HLGR using canned message (See Blood and Sterile Fluids HLGR Results Reporting).

⁴ Report only if Pen is I or R

⁵ Report on PMH patients; base on Ampicillin result

⁶ E. gallinarum and E. casseliflavus report as \mathbf{R} with the statement: "This organism always has intrinsic non-transmissible resistance to vancomycin. The patient does not require isolation. ⁷ If Vancomycin and Ampicillin are R except *E. gallinarum* and *E. casseliflavus*.

⁸ Base on E-test

^a includes MSH, PMH, Baycrest, TRI, CAMH

^b includes TG, TW, Bridgepoint, Grace

^c includes CHC, Ajax/Pickering

Note: Listeria species - DO NOT report susceptibility result. Report with ISOLATE comment - "Routine in vitro susceptibility testing of this organism is unreliable. *Listeria* spp. should be considered resistant to all cephalosporins. The recommended regimen for therapy is ampicillin. If additional advice on antimicrobial therapy is required, please contact the Medical Microbiologist."

Corynebacterim species, Bacillus species - DO NOT report susceptibility result. Report with ISOLATE comment "In vitro susceptibility testing for this organism is not routinely performed and/or is unreliable. If advice on antimicrobial therapy is required, please contact the Medical Microbiologist".

If all antimicrobial agents are resistant, inform the Microbiologist on-call

TML\MSH Microbiology Department	Policy #MI\ANTI\03\10\v10	Page 1 of 1			
Policy & Procedure Manual					
Section: Antimicrobial Susceptibility Testing	Subject Title: Spinal Fluids - G	ram Negative			
Manual	Susceptibility Reporting - 1				
Issued by I ADODATODY MANACED	Original Datas January 10, 2000				
Issued by: LADOKATOKI MANAGEK	Original Date: January 10, 2000				
Approved by: Laboratory Director	Revision Date: November 21, 20	005			

Spinal Fluids – Gram Negative Susceptibility Reporting – 1 – Enterobacteriaceae and Acinetobacter spp., Salmonella species including S. typhi

Antimicrobial Agent	Enterobacter	iaceae and Acin	etobacter spp.	Salmonella species		
		k k	including S. typhi			
	MSH ^a	UHN ^b	CHC	MSH ^a	UHN⁰	CHC
Amikacin	X′	X'	X′			
Ampicillin	$X^{4, 6}$	X ^{4, 6}	X ^{4, 6}	Х	Х	Х
Ceftazidime	$X^{2, 6}$	$X^{2, 6}$	$X^{2, 6}$			
Ceftriaxone	X ^{5, 6}	X ^{5, 6}	X ^{5, 6}	Х	Х	Х
Colistin	X ⁸	X ⁸	X ⁸			
Gentamicin	X	Х	Х			
Imipenem	X ³	X^3	X^3			
Piperacillin/Tazobactam	X°	X ⁵				
Trimethoprim/Sulfa	X	X	X	Х	X	X
Tobramycin	X	X	X			

¹ Report if Genta R

^{2} Report only if R.

³ Report if R **OR** if R to **All** other Antimicrobial Agents **OR** if only aminoglycoside is S

- ⁴ Always report Klebsiella spp., Enterobacter spp., Citrobacter spp., Proteus vulgaris, Proteus penneri, Acinetobacter spp., H. alvei and Serratia species as Ampicillin R
- ⁵ Cedecea spp., Citrobacter spp., Enterobacter spp., Hafnia spp., Morganella morganii, Pantoea agglonerans, Proteus penneri, Proteus vulgaris, Providencia species, Seratia species, if S, report with comment "Resistance to extended-spectrum penicillins, beta-lactam/beta-lactamase inhibitor combinations, and cephalosporins may develop during therapy with these agents. For serious infections, these agents should be avoided and consultation with a medical microbiologist or infectious disease physician is strongly recommended."

⁶ For *Enterobacteriaceae* if any one of cefotaxime/ceftriaxone or ceftazidime=R, report all as R

⁷ Report if both Gentamicin and Tobramycin are R.
 ⁸ For Acinetobacter species, report if R to all other drugs; base on Polymyxin B.

^a includes MSH, PMH, Baycrest, TRI, CAMH

^b includes TG, TW, Bridgepoint, Grace

^c includes CHC, Ajax/Pickering

Note: If all antimicrobial agents are resistant, inform the Microbiologist on-call

Page 20

TML\MSH Microbiology Department Policy & Procedure Manual	Policy #MI\ANTI\03\11\v10	Page 1 of 1			
Section: Antimicrobial Susceptibility Testing	Subject Title: Spinal Fluids - Gram Negative				
Manual	Susceptibility Reporting – 2				
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000				
Approved by: Laboratory Director	Revision Date: November 21, 20)05			
	Review Date: April 21 2005				

Spinal Fluids – Gram Negative Susceptibility Reporting – 2 – *Pseudomonas aeruginosa*, *Pseudomonas* spp. (other than P. aeruginosa), fastidious gram-negative bacteria, nonfermenters, M. catarrhalis, N. gonorrhoeae, N. meningitidis, Stenotrophomonas maltophilia, Burkholderia cepacia, Haemophilus species.

Antimicrobial Agent	Ps. aeruginosa		S. maltophilia	B. cepacia	Haemophilus species.	
	MSH ^a	UHN ^b	CHC ^c	All sites	All sites	All sites
Amikacin	X^3	X	X^3			
Ampicillin	R					X ^I
Ceftazidime	Х	Х	Х		Х	
Ceftriaxone						Х
Colistin	X^4	X^4	X^4	X ⁴	X ⁴	
Gentamicin	Х	Х	Х			
Imipenem	X^2	X^2	X^2			
Meropenem					Х	
Piperacillin/Tazobactam	Х	Х	Х			
TMP/SMX	R			Х	Х	
Tobramycin	X	X	X			

¹ Base on beta-Lactamase result and KB Ampicillin

² Report if R **OR** if R to **All** other Antimicrobial Agents **OR** if only aminoglycoside is S

³Report if both Gentamicin and Tobramycin are R.

⁴Report if R to all other drugs; base on Polymyxin B.

^a includes MSH, PMH, Baycrest, TRI, CAMH

^b includes TG, TW, Bridgepoint, Grace

^c includes CHC, Ajax/Pickering

Note: *Pseudomonas* species (other than *P. aeruginosa*), fastidious gram-negative bacteria, nonfermenters, *N. gonorrhoeae* and *N. meningitidis* - DO NOT report susceptibility result. Report with ISOLATE comment "In vitro susceptibility testing for this organism is not routinely performed and/or is unreliable. If advice on antimicrobial therapy is required, please contact the Medical Microbiologist".

For *M. catarrhalis* - DO NOT report susceptibility result. Report with ISOLATE comment: "The majority of *Moraxella catarrhalis* are resistant to ampicillin. In vitro susceptibility testing for this organism is not routinely performed and/or is unreliable. If advice on antimicrobial therapy is required, please contact the Medical Microbiologist".

If all antimicrobial agents are resistant, inform the Microbiologist on-call.

TML\MSH Microbiology Department Policy & Procedure Manual	Policy #MI\ANTI\03\12\v10	Page 1 of 1		
Section: Antimicrobial Susceptibility Testing	g Subject Title: Blood and other Sterile Sites -			
Manual	Gram Positive Susceptibility			
	Reporting – 1			
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000			
Approved by: Laboratory Director	Revision Date: November 21, 20)05		
	Review Date: April 21, 2005			

Blood and other Sterile Sites - Gram Positive Susceptibility Reporting – 1 – *Staphylococcus, Enterococcus*

Antimicrobial Agent	Staph	ylococcus s	pecies	Enterococcus species			
	MSH ^a	UHN [⊳]	CHC ^c	MSH ^a	UHN⁵	CHC ^c	
Ampicillin				Х	Х	Х	
Cefazolin	X^2	X^2	X^2				
Clindamycin	X ³	X ⁵	X°				
Cloxacillin	X^2	X^2	X^2				
Erythromycin	X ⁵	X ⁵	X ³				
HLGR ³				Х	Х	Х	
HLSR ³				X^4	X^4	X ⁴	
Linezolid				X′	Χ'	Χ'	
Pipercillin/Tazobactam				X ¹	X ⁱ		
Synercid				X′	Χ'	Χ′	
TMP/SMX	Х	Х	Х				
Vancomycin	X ⁸	Х	Х	X ⁶	X ⁶	X ⁶	

¹ Base on ampicillin result (PMH, TGH, TWH)

² Base on Oxacillin result

 ³ HLGR = High Level Gentamicin Resistant; HLSR = High Level Streptomycin Resistant Report base on HLGR using canned message (See Blood and Sterile Fluids HLGR Results Reporting).

⁴ Report only if requested.

⁵ Clinda R = Erythro R

 6 *E. gallinarum and E. casseliflavus* report as **R** with the statement "This organism always has intrinsic non-transmissible resistance to vancomycin. The patient does not require isolation."

⁷ If Vancomycin and Ampicillin are R except *E. gallinarum* and *E. casseliflavus*.

⁸ Only if Oxacillin=R.

^a includes MSH, PMH, Baycrest, TRI, CAMH

^b includes TG, TW, Bridgepoint, Grace

^c includes CHC, Ajax/Pickering

Note: If all antimicrobial agents are resistant, inform the Microbiologist on-call.

TML\MSH Microbiology Department	Policy #MI\ANTI\03\13\v10	Page 1 of 1			
Policy & Procedure Manual					
Section: Antimicrobial Susceptibility Testing	Subject Title: Blood and other Sterile Sites -				
Manual	Gram Positive Susceptibility				
	Reporting – 2				
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000				
Approved by: Laboratory Director	Revision Date: November 21, 20)05			
	Review Date: April 21, 2005				

Blood and other Sterile Sites - Gram Positive Susceptibility Reporting – 2 - S. pneumoniae, viridans Streptococcus, Streptococcus bovis,, S. milleri group, Group A, B, C, G Streptococcus, Listeria species, Corynebacterim species, Bacillus species

Antimicrobial	<i>S. p</i>	neumor	niae	viridans <i>Strep</i> .			S. milleri	Group A, B, C,G
Agent			Streptococcus bovis		group	Strep.		
	MSH ^a	UHN ^b	CHC ^c	MSH ^a	UHN [♭]	CHC ^c	All Sites	All Sites
Ceftriaxone	X ^{2, 5}	X ^{2, 5}	X ^{2, 5}	X ^{2, 5}	X ^{2, 5}	X ^{2, 5}	X ^{2,5}	
Cefepime	X ⁶			X ⁶				
Clindamycin	Х	Х	Х					X ³
Erythromycin	Х	Х	Х					X ³
Levofloxacin	X ^{1, 8}		X ¹					
Moxifloxacin	X ^{1,7}	X ^{1,7}						
Penicillin	X^2	X^2	X^2	$X^{2,4}$	$X^{2, 4}$	X ^{2, 4}	$X^{2,4}$	Х
Vancomycin	X ⁵	X ⁵	X ⁵	X ⁵	X^5	X ⁵	X ^{2, 5}	X ⁵

¹Adults only (>18 y)

² Base on E-test

³ Report as R if D-zone is present. Report clindamycin as R if erythromycin is R

⁵ Report only if Pen I or R

⁴ For viridans Streptococcus and S. milleri, also report MIC value as Isolate Comment

⁶ Report on PMH patients only if Ceftriaxone is I or R, base on Ceftriaxone result.

⁷ Base on Levofloxacin result; Report on MSH, PMH, TG and TW patients only.

⁸ DO NOT report on MSH and PMH patients.

^a includes MSH, PMH, Baycrest, TRI, CAMH

^b includes TG, TW, Bridgepoint, Grace

^c includes CHC, Ajax/Pickering

Note: *Listeria* species – DO NOT report susceptibility result. Report with ISOLATE comment – "Routine in vitro susceptibility testing of this organism is unreliable. *Listeria* spp. should be considered resistant to all cephalosporins. The recommended regimen for therapy is ampicillin. If additional advice on antimicrobial therapy is required, please contact the Medical Microbiologist."

Corynebacterim species, *Bacillus* species. - DO NOT report susceptibility result. Report with ISOLATE comment "In vitro susceptibility testing for this organism is not routinely performed and/or is unreliable. If advice on antimicrobial therapy is required, please contact the Medical Microbiologist".

If all antimicrobial agents are resistant, inform the Microbiologist on-call.

TML\MSH Microbiology Department	Policy #MI\ANTI\03\14\v10	Page 1 of 1			
Policy & Procedure Manual					
Section: Antimicrobial Susceptibility Testing	Antimicrobial Susceptibility Testing Subject Title: Blood and other Sterile Sites				
Manual	Gram Negative Susceptibility				
		= -			
	Reporting – 1				
Issued by: LABORATORY MANAGER	Reporting – 1 Original Date: January 10, 2000				
Issued by: LABORATORY MANAGER Approved by: Laboratory Director	Reporting – 1 Original Date: January 10, 2000 Revision Date: November 21, 20	005			

Blood and other Sterile Sites - Gram Negative Susceptibility Reporting -1 - Enterobacteriaceae and *Acinetobacter* spp., *Salmonella* species including *S. typhi*

	Enterobacteriaceae and Acinetobacter spp.			Salmonella spp. including S. typ		
Antimicrobial Agent	MSH ^a	UHN ^b	CHC ^c	MSH ^a	UHN ^b	CHC
Amikacin	X^{13}	X ¹³	X^{13}			
Ampicillin	X ¹	X ¹	X ¹	Х	Х	Х
Cefazolin	$X^{2, 12}$	$X^{2, 12}$	X^2			
Cefepime	X ^{5, 12}					
Ceftazidime	X ^{6, 9, 12}	X ^{6, 12}	$X^{6, 12}$			
Ceftriaxone	$X^{10, 12}$	$X^{10, 12}$	$X^{10, 12}$	Х	Х	Х
Ciprofloxacin	X ⁴	X ⁴	X ⁴	X4	X^4	X^4
Colistin	X ¹⁴	X ¹⁴	X^{14}			
Gentamicin	X	Х	Х			
Imipenem	X ⁸	X ⁸	X ⁸			
Nalidixic Acid				X ¹¹	X^{II}	X^{II}
Piperacillin/Tazobactam	X ¹⁰	X ¹⁰				
TMP/SMX	X	X	X	Х	Х	Х
Tobramycin	X	$X^{3,7}$	X'			

¹ Always report Klebsiella spp, Enterobacter spp, Acinetobacter spp., H. alvei, Citrobacter spp., Proteus vulgaris, Proteus penneri, and Serratia species as R

² Always report Enterobacter spp, Citrobacter spp., Proteus vulgaris, Proteus penneri, Acinetobacter spp, and Serratia species as R.

³ TG and TW – report Tobramycin on all PD effluents

⁴ Adults only (>1 $\overline{8}$ y)

⁵ Report on PMH patients only if Cetriaxone is I or R.

⁶ Report only if R

⁷ Report if Genta is R

⁸ Report if R **OR** if R to **All** other Antimicrobial Agents **OR** if only aminoglycoside is S

⁹ Always report for PMH

¹⁰ Cedecea spp., Citrobacter spp., Enterobacter spp., Hafnia spp., Morganella morganii, Pantoea agglonerans, Proteus penneri, Proteus vulgaris, Providencia species, Seratia species, if **S**, report with comment "Resistance to extended-spectrum penicillins, beta-lactam/beta-lactamase inhibitor combinations, and cephalosporins may develop during therapy with these agents. For serious infections, these agents should be avoided and consultation with a medical microbiologist or infectious disease physician is strongly recommended."

¹¹DO NOT Report. If Ciprofloxacin is S and Nalidixic Acid is R, report isolate comment "This isolate is tested Fluoroquinolone-susceptible and nalidixic acid-resistant; may be associated with clinical failure or delayed response in fluoroquinolone-treated patients with extraintesinal salmonellosis."

¹² For *Enterobacteria ceae* if any one of cefotaxime/ceftriaxone or ceftazidime=R, report all as R

¹³Report if both Gentamicin and Tobramycin are R.

¹⁴ For Acinetobacter species, report if R to all other drugs; base on Polymyxin B.

^a includes MSH, PMH, Baycrest, TRI, CAMH

^b includes TG, TW, Bridgepoint, Grace

^c includes CHC, Ajax/Pickering

Note: If all antimicrobial agents are resistant, inform the Microbiologist on-call

PROCEDURE MANUAL

TML\MSH Microbiology Department	Policy #MI\ANTI\03\15\v10	Page 1 of 1			
Policy & Procedure Manual	-				
Section: Antimicrobial Susceptibility Testing	Subject Title: Blood and other Sterile Sites - Gram				
Manual	Negative Susceptibility Reporting – 2				
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000				
Approved by: Laboratory Director	Revision Date: November 21, 2005				
	Review Date: April 21, 2005				

Blood and other Sterile Sites - Gram Negative Susceptibility Reporting -2 - *Pseudomonas* aeruginosa, *Pseudomonas* spp. (other than P. aeruginosa), fastidious gram-negative bacteria, non-fermenters, *M. catarrhalis*, *N. gonorrhoeae*, *N. meningitidis*, *Stenotrophomonas maltophilia*, *Burkholderia cepacia*, *Haemophilus* species.

Antimicrobial Agent	P. aeruginosa			S. maltophilia	B. cepacia	Haemophilus species
	MSH ^a	UHN⁵	CHC ^c	All sites	All sites	All sites
Amikacin	X ⁵	X ⁵	X ⁵			
Ampicillin	R					X
Cefazolin	R					
Cefepime	X^3					
Ceftazidime	Х	Х	Х		Х	
Ceftriaxone						Х
Ciprofloxacin	X^2	X^2	X^2			\mathbf{X}^2
Colistin	X^6	X^6	X^{6}	X^{6}	X ⁶	
Gentamicin	Х	Х	Х			
Imipenem	X^4	X^4	X^4			
Levofloxacin				X^2		
Meropenem					Х	
Piperacillin/Tazobactam	Х	Х	Х			
TMP/SMX	R			Х	Х	
Tobramycin	Х	Х	Х			

¹ Base on beta -lactamase result and KB Ampicillin

² Adults only (>18 y)

³ For PMH only

⁴ Report if R **OR** if R to **All** other Antimicrobial Agents **OR** if only aminoglycoside is S

⁵ Report if both Gentamicin and Tobramycin are R.

⁶ Report if R to all other drugs; base on Polymyxin B.

^a includes MSH, PMH, Baycrest, TRI, CAMH

^b includes TG, TW, Bridgepoint, Grace

^c includes CHC, Ajax/Pickering

Note: Pseudomonas species (other than P. aeruginosa), fastidious gram-negative bacteria, non-

fermenters, *N. gonorrhoeae* and *N. meningitidis* - DO NOT report susceptibility result. Report with ISOLATE comment "In vitro susceptibility testing for this organism is not routinely performed and/or is unreliable. If advice on antimicrobial therapy is required, please contact the Medical Microbiologist".

For *M. catarrhalis* - DO NOT report susceptibility result. Report with ISOLATE comment: "The majority of *Moraxella catarrhalis* are resistant to ampicillin. In vitro susceptibility testing for this organism is not routinely performed and/or is unreliable. If advice on antimicrobial therapy is required, please contact the Medical Microbiologist".

If all antimicrobial agents are resistant, inform the Microbiologist on-call

TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI\ANTI\04\v01	Page 1 of 5	
Section: Antimicrobial Susceptibility Testing	Subject Title: Antimicrobial Related LIS		
Manual	Canned Messages		
Issued by: LABORATORY MANAGER	Original Date: November 21, 2005		
Approved by: Laboratory Director	Revision Date:		

ANTIMICROBIAL RELATED LIS CANNED MESSAGES

I. <u>Introduction</u>

Antimicrobial related canned messages are built into the Laboratory Information System to provide uniform reporting phrases to be used when certain pre-described conditions are met.

II. <u>Procedure</u>

A. Automatic Canned Messages:

The lists below are automatic canned messages that will appear when set conditions are met. The message will appear in a warning box before exiting an order.

- 1. When the message code appears, press F12 to save.
- 2. Continue with another F12 to save the order.
- 3. View the report.
- 4. If the same message had been save previously (i.e. appeared more than once), go to the Isolate Comment window, press CTLR and L at the line to remove the duplicate line.
- 5. Re-status as required.
- 6. Press F12 to save the order.

Ear and Eye specimens with susceptibility results

LIS Isolate Canned Message Code: &eye; attached to Organism classes A and B with sources EYE, EYEN, EAR and drugs am, betalac, cc, peng, sxt.

"These susceptibility testing results are based on guidelines for systemic antimicrobial agents and may not accurately represent activity of topical agents."

For MSH MRO's

LIS Isolate Canned Message Code: \MRES, attached to drug – tax "MULTIPLY ANTIBIOTIC RESISTANT ORGANISM. THIS PATIENT MUST BE ON "CONTACT PRECAUTIONS" UNTIL FURTHER NOTICE FROM INFECTION CONTROL."

For MSH MRSA's

LIS Isolate Canned Message Code: \MRSI, attached to organism staamr "THIS PATIENT MUST BE ON "MRSA PRECAUTION" UNTIL FURTHER NOTICE."

MRSA isolated from MRSA Screen Susceptibility Result Comment

LIS Isolate Canned Message Code: **MRSS**; linked to drug code - oxa "Susceptibility results are provided to guide decolonization therapy; if decolonization therapy is being considered, please consult your infection control team."

TML\MSH Microbiology Department	Policy # MI\ANTI\04\v01	Page 2 of 5
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual	Antimicrobial Related LIS Canned	
	Messages	

Staphylococcus saprophyticus and CNST from urine

LIS Isolate Canned Message Code: \ssap; attached to organism stasap. "Susceptibility testing of this organism is not warranted because infections respond to concentrations achieved in urine of antimicrobial agents commonly used to treat acute, uncomplicated urinary tract infections e.g. nitrofurantoin, trimethoprim/sulfa or fluoroquinolones."

β-haemolytic Streptococcus Groups A, B, C and G

LIS Isolate Canned Message Code: \GBS; attached to organism straga, strpyo, strgrc, strgrg.

"This organism is intrinsically susceptible to penicillin. If treatment is required AND this patient cannot be treated with penicillin, please contact the Microbiology Department within 48 hours to request sensitivity testing."

For MSH VRE's

LIS Isolate Canned Message Code: \VREI, attached to organisms entfac and entfae "THIS PATIENT MUST BE ON "VRE PRECAUTION" UNTIL FURTHER NOTICE."

Vancomycin for *E. gallinarum, and E. casseliflavus* report as **R** with the statement LIS Isolate Canned Message Code: **\EntV**; attached to organisms - entgal and entcas. "This organism always has intrinsic non-transmissible resistance to vancomycin. The patient does not require isolation."

For Listeria species:

LIS Isolate Canned Message Code: \LIST; attached to organisms lismoc and lismon. "Routine in vitro susceptibility testing of this organism is unreliable. *Listeria* spp. should be considered resistant to all cephalosporins. The recommended regimen for therapy is ampicillin. If additional advice on antimicrobial therapy is required, please contact the Medical Microbiologist."

For isolates that for which **susceptibility testing is not routinely performed and/or is unreliable**:

LIS Isolate Canned Message Code: \NSEN; attached to organisms and Isolate Comment keypad.

"In vitro susceptibility testing for this organism is not routinely performed and/or is unreliable. If advice on antimicrobial therapy is required, please contact the Medical Microbiologist".

TML\MSH Microbiology Department	Policy # MI\ANTI\04\v01	Page 3 of 5
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual	: Antimicrobial Related LIS Canned	
	Messages	

For *Cedecea* spp., *Citrobacter* spp., *Enterobacter* spp., *Hafnia* spp., *Morganella morganii*, *Pantoea agglonerans*, *Proteus penneri*, *Proteus vulgaris*, *Providencia* species, *Serratia* species that are ceftriaxone or pipercillin/tazobactam = S: LIS Isolate Canned Message Code: &cit or &ent or &ser, attached to Organisms ceddav, cedlap, cedspp, Classes d, e, H, L, S, and T.

"Resistance to extended-spectrum penicillins, beta-lactam/beta-lactamase inhibitor combinations (e.g. pipercillin/tazobactam), and cephalosporins may develop during therapy with these agents. For serious infections, these agents should be avoided and consultation with a medical microbiologist or infectious disease physician is strongly recommended."

For **MSH** *E. coli, Klebsiella* species, *Proteus* Class A ESBL, Infection Control message: Isolate canned message code **&taz** linked to organisms *E. coli*, Class J and Class T: "MULTIPLY ANTIBIOTIC RESISTANT ORGANISM. THIS PATIENT MUST BE ON "CONTACT PRECAUTIONS" UNTIL FURTHER NOTICE FROM INFECTION CONTROL."

ESBL Comments

Attached to organisms esccol, Class J and Class T

Desbinh=Y Dfox=S **~&cla** "The susceptibility pattern suggests that this organism contains a class A extended spectrum beta-lactamase (ESBL)."

Dtaz=R Desbinh=N Dfox=R or I ~&claC "The susceptibility pattern suggests that this organism contains a class C extended spectrum beta-lactamase (ESBL)."

Desbinh=Y Dfox=R or I **~&clAC** 'The susceptibility pattern suggests that this organism contains class A and C extended spectrum beta-lactamases (ESBL)."

Dtaz=R Desbinh=N Dfox=S **~&esbl** "The susceptibility pattern suggests that this organism contains an extended spectrum beta-lactamase (ESBL) other than class A or C."

For *Haemophilus* species from Respiratory and Miscellaneous Sites – LIS Isolate Canned Message Code, attached to organism Class X:

If beta-lactamase is negative, **BLa-** "beta-lactamase negative result suggests susceptible to ampicillin."

If beta-lactamase is positive, **BLa**+ "beta-lactamase positive result suggests resistance to ampicillin but generally susceptible to amoxicillin-clavulanic acid and cefuroxime."

TML\MSH Microbiology Department	Policy # MI\ANTI\04\v01	Page 4 of 5
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual	Antimicrobial Related LIS Canned	
	Messages	

For *Salmonella* species:

Nalidixic acid confirming Fluoroquinolone Susceptibility in *Salmonella* species Preliminary Report Comment

"Susceptibility testing suggests that this organism may not be fully susceptible to fluoroquinolones. Further testing is being completed and results will be available tomorrow. If you are interested in preliminary susceptibility results, please contact the microbiology laboratory."

Comment after KB – Nalidixic Acid has been completed:

LIS Isolate Canned Message Code: **&sal**; attached to organisms Class q and saltyp: Dcip=S Dna=R RI RR R1 Ra Rb Re Rl Rk ~&sal

If Ciprofloxacin is S and Nalidixic Acid is R, report isolate comment "This isolate is tested Fluoroquinolone-susceptible and nalidixic acid-resistant; may be associated with clinical failure or delayed response in fluoroquinolone-treated patients with extraintesinal salmonellosis."

For *M. catarrhalis* - LIS Isolate Canned Message Code: \mcat; attached to Organism: "The majority of *Moraxella catarrhalis* are resistant to ampicillin. In vitro susceptibility testing for this organism is not routinely performed and/or is unreliable. If advice on antimicrobial therapy is required, please contact the Medical Microbiologist".

B. Canned Messages to be selected from the Isolate Comment keypad:

The listed below are canned messages to be selected from the Isolate Comment keypad when needed.

- 1. At the LIS Isolate Comment Window, type the appropriate number on the keypad.
- 2. Press F12 to save.
- 3. Continue with another F12 to save the order.
- 4. View the report.
- 5. Status the report as required.

For **INH** reporting if 0.1 mg/L=R and 0.4mg/L=S:

LIS Isolate Canned Message Code: **INHr**; select from keypad "This isolate has low-level resistance to isoniazid (INH). Patients infected with strains exhibiting this level of INH resistance may benefit from continuing therapy with INH. Consultation with a specialist experienced in the treatment of tuberculosis is recommended."

BORSA (DENKA *mec* A-negative *S. aureus* with oxacillin MIC>=4mg/L)

LIS Isolate Canned Message Code: **BORS**; select from Isolate Keypad "This organism is resistant to cloxacillin by a mechanism different from that in typical MRSA. Consultation with a Microbiologist or Infectious Disease physician is advised."

TML\MSH Microbiology Department	Policy # MI\ANTI\04\v01	Page 5 of 5
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual	: Antimicrobial Related LIS Canned	
	Messages	

If susceptibility is done on request for *β-haemolytic Streptococcus* Groups A, B, C & G Do not remove original canned message. Add message from Isolate Comment keypad code "Susceptibility\done" – "Susceptibility completed as requested."

Enterococcus from Blood and Sterile Sites:

If high level gentamicin is **susceptible** (regardless of streptomycin result), select from Isolate Comment Keypad **\EGMS**: "Serious enterococcal infections may require an aminoglycoside for synergy. Please contact the Medical Microbiologist for treatment advice".

If high level gentamicin is **resistant** (regardless of streptomycin result), select from Isolate Comment Keypad **\EGMR**: "This organism is high level aminoglycoside resistant. Please contact the Medical Microbiologist for treatment advice".

E. coli, Klebsiella species, *Proteus* cefpodoxime=I or **R**, set up KB ESBL, Preliminary Report Comment; select from ISOLATE keypad - \podR

"Susceptibility testing suggests that this organism may contain an extended spectrum betalactamase. Further testing is being completed and results will be available tomorrow. If you are interested in preliminary susceptibility results, please contact the microbiology laboratory."

TML\MSH Microbiology Department	Policy #MI\ANTI\04\01\v03	Page 1 of 2
Policy & Procedure Manual		
Section: Antimicrobial Susceptibility Testing	Subject Title: Appendix I – Disk Diffusion	
Manual		
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000	
Approved by: Laboratory Director	Revision Date: November 21, 20)05

APPENDIX I - DISK DIFFUSION

I. Introduction

The disk diffusion method of susceptibility testing (also known as the Kirby-Bauer (KB) method) has been standardized primarily for testing of rapidly growing bacteria. To perform the test, filter paper disks impregnated with a specific amount of antimicrobial agent are applied to the surface of an agar medium that has been inoculated with a known amount of the test organism. The drug in the disk diffuses through the agar. As the distance from the disk increases, the concentration of the antimicrobial agent decreases creating a gradient of drug concentrations in the agar medium. Concomitant with diffusion of the drug, the bacteria that were inoculated and that are not inhibited by the concentration of the antimicrobial agent continue to multiply until a lawn of growth is visible. In areas where the concentration of drug is inhibitory, no growth occurs, forming a zone of inhibition around each disk. Criteria currently recommended for interpreting zone diameters and MIC results for commonly used antimicrobial agents are published by NCCLS. Results are reported categorically as Susceptible (S), Intermediate (I), or Resistant (R).

II. <u>Materials</u>

Antimicrobial disks (store frozen with a desiccant) Mueller Hinton Agar (MH) Mueller Hinton Blood Agar (MHB) Haemophilus Test Media (HTM) Trypticase Soy Broth (TSB) (3 mL) VITEK colorimeter Sterile saline Sterile swabs

III. Procedure

- 1. Allow disks to come to room temperature before opening the container.
- 2. Using the Vitek colorimeter, prepare a suspension of the test organism in sterile saline equivalent to a 0.5 McFarland standard using isolated colonies. If there is not enough growth, inoculate the organism into TSB, and incubate at 35°C for 2-4 hours or until it reaches the turbidity of a 0.5 McFarland standard.

TML\MSH Microbiology Department	Policy # MI\ANTI\04\01\v03	Page 2 of 2
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual		

- 3. Using a sterile cotton swab, inoculate the organism onto an appropriate agar plate, streaking in 3 directions over the entire agar surface. For organisms that grow rapidly use MH agar. For *Haemophilus* species use HTM and for *S. pneumoniae* use MHB. For other organisms that do not grow on MH, use MHB.
- 4. Using forceps or a disk dispenser, apply the appropriate Antimicrobial disks onto the agar. Place the disks with an equal distance apart from each other and put no more than 6 disks on a 100mm diameter plate.
- 5. Incubate plates as follows: *Campylobacter* species - microaerophilically at 35°C x 18 hours *Haemophilus* species - CO₂, 35°C x 18 hours *S. pneumoniae* - CO₂, 35°C x 20 to 24 hours *S. aureus* and *Enterococcus* species for Methicillin and Vancomycin - O₂, 35°C x 24 hours Others - O₂, 35°C x 18 hours

IV. Interpretation

After incubation, measure the diameters of the zone of complete inhibition with callipers using transmitted light. For MH agar, measure from the back of the plate.

Refer to NCCLS Document M100-S15 (M2) for the zone size interpretations. Report susceptible, resistant and intermediate as appropriate.

V. Quality Control

Test the following organisms each time a new batch of MH agar is prepared and once weekly. Subculture the organisms from the TSB slant (in fridge) to BA the day before setting up the QC. *S. aureus* ATCC 25923

E. coli ATCC 25922 *P. aeruginosa* ATCC 27853 In addition, test *S. faecalis* ATCC 29212 each time a new batch of MH is prepared.

Perform weekly Quality Control on HTM agar with *Haemophilus influenzae* ATCC 49247. Test for growth of *Haemophilus influenzae* ATCC 10211 on each new batch HTM.

See NCCLS Document M100-S15 (M2) Table 3 for acceptable QC results.

VI. <u>Reference</u>

Clinical and Laboratory Standards Institute (CLSI) Document - Performance Standards for Antimicrobial Disk Susceptibility Testing M2-A8, 2003.

TML\MSH Microbiology Department Policy & Procedure Manual	Policy #MI\ANTI\04\02\v01	Page 1 of 4
Section: Antimicrobial Susceptibility Testing	Subject Title: Appendix II – Double Disk	
Manual	Diffusion for Erythromycin and	
	Clindamycin on ß-haemolytic Streptococci	
	and Streptococcus pneumoniae	
Issued by: LABORATORY MANAGER	Original Date: November 21, 20	05
Approved by: Laboratory Director	Revision Date:	

APPENDIX II – DOUBLE DISK DIFFUSION for ERYTHROMYCIN and CLINDAMYCIN on & haemolytic Streptococci Groups A, B, C, G and Streptococcus pneumoniae

I. <u>Introduction</u>

Macrolide (erythromycin) resistant ß-haemolytic *Streptococci* and *Streptococcus pneumoniae* isolates may have constitutive or inducible resistance to lincosamides (clindamycin). The mechanisms of resistance include:

- Ribosomal modification encoded by an *erm* gene; also refer to as MLS_B (macrolide, lincosamide and type B streptogramin) resistance.
- Efflux of the antibiotic encoded by a *mef* gene; resistant only to macrolide
- Drug inactivation

Inducible clindamycin resistance can be detected using a disk approximation test with a clindamycin disk placed beside an erythromycin disk as part of the normal disk diffusion test.

II. <u>Materials</u>

Antimicrobial disks – clindamycin (DA, 2 μ g) and erythromycin (E, 15 μ g) Mueller Hinton Blood Agar (MHB) Trypticase Soy Broth (TSB) (3 mL) VITEK colorimeter Sterile saline Sterile swabs

III. <u>Procedure</u>

- 1. Allow disks to come to room temperature before opening the container.
- 2. Using the Vitek colorimeter, prepare a suspension of the test organism in sterile saline equivalent to a 0.5 McFarland standard using isolated colonies.
- 3. Using a sterile cotton swab, inoculate the standardized organism onto an MHB agar plate, streak in three directions over the entire agar surface.
- 4. Place plate on disk template (Figure 1.)
- 5. Using forceps or a disk dispenser, apply the clindamycin and erythromycin disks onto the agar 12 mm apart from each other edge to edge using template below (Figure 1). Other antimicrobial disks can be placed on the same agar plate if needed.



Figure 1. Template for Clindamycin and Erythromycin disks placement

6. Incubate plates in CO_2 at $35^{\circ}C$ for 20 to 24 hours

IV. <u>Interpretation</u>

- After incubation, measure the diameters of the zone of complete inhibition with callipers/ruler. Refer to Clinical and Laboratory Standards Institute (CLSI) Document -M100-S15 (M2) for the zone size interpretations.
- 2. Enter zone size measurements into the LIS.
- 3. Organisms that show flattening of the clindamycin zone adjacent to the erythromycin disk in the shape of the letter D (referred to as a "D" zone) have inducible clindamycin resistance. Enter into the LIS the presence or absence of "D" zone as "Y" or "N" under LIS drug "D zone". Isolates that show the presence of D zone will be automatically reflexed in the LIS to report as "clindamycin resistant".

TML\MSH Microbiology Department	Policy # MI\ANTI\04\02\v01	Page 3 of 4
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual		

Examples of Zone of Inhibition Patterns and their Interpretation





Page 35

TML\MSH Microbiology Department	Policy # MI\ANTI\04\02\v01	Page 4 of 4
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual		

V. <u>Quality Control</u>

See Clinical and Laboratory Standards Institute (CLSI) Document - M100-S15 (M2) Table 3 for acceptable QC results.

VI. <u>References</u>

Clinical and Laboratory Standards Institute (CLSI) Document - Performance Standards for Antimicrobial Disk Susceptibility Testing M2-A8, 2003.

Clinical and Laboratory Standards Institute (CLSI) Document - Performance Standards for Antimicrobial Disk Susceptibility Testing Information Supplement Table 2H M2-Disk Diffusion M100-S15 2005.

Quality Management Program-Laboratory Services (QMP-LS) Committee Comments BACT-020, Vol. 3, 2.2:721-724.

Streptococci and *Staphylococcus* (overview of macrolides and lincosamide resistance) Leclercq. CID 2002;34:482-92

Streptococcus pneumoniae Descheemaeker et al. JAC 2000 45:167-173

Beta-haemolytic streptococcus (Groups A,B,C,G) GAS Descheemaeker et al. JAC 2000 45:167-173 GBS de Azavedo et al. AAC 1001;45:3504-3508 GCS & GGS Kataja et al. AAC 1998;42:1493-1494

TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI\ANTI\04\03\v03	Page 1 of 5
Section: Antimicrobial Susceptibility Testing	Subject Title: Appendix III - Do	ouble Disk
Manual	Test for ESBL	
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000	
Approved by: Laboratory Director	Revision Date: November 21, 20)05

APPENDIX III - DOUBLE DISK TEST FOR ESBL

I. Introduction

Class A or Bush Group 1 extended spectrum beta-lactamases (ESBLs) are inhibited by clavulanic acid. This may be detected by testing the suspected organism to a 3rd generation cephalosporin alone and in combination with clavulanic acid. If the combination results in an expanded zone of inhibition compared to that of the 3rd generation cephalosporin alone, it is indicative of the presence of an ESBL.

II. <u>Materials</u>

Mueller-Hinton (MH) agar (150) mm 20/10 mg amoxicillin-clavulanate disc 30 mg ceftazidime disc 30 mg ceftriaxone or cefotaxime disc 30 mg aztreonam disc 10 mg cefpodoxime disc (optional) 30 mg cefoxitin disc Quality control strain: *E. coli* ATCC 51446

III. Procedure

- 1. Prepare a bacterial suspension of the organism to be tested that has a turbidity equivalent to that of a 0.5 McFarland standard.
- 2. Inoculate a Mueller-Hinton agar plate with this suspension in accordance with NCCLS M100-S10 (M2) guidelines for disc diffusion testing.
- 3. Place the amoxicillin-clavulanic acid disc towards the centre of the plate.
- 4. Carefully measure 15 mm out from the edge of that disc at 90° angles marking the plate.
- 5. Place a ceftazidime disk on the plate so that its inner edge is 15 mm (the mark) from the amoxicillin-clavulanic acid disc (See Figure 1 KB-ESBL Template).
- 6. Do the same with cefotaxime (or ceftriaxone), aztreonam and cefpodoxime discs so that they are spaced 90° apart and 15mm from the centre disc.

TML\MSH Microbiology Department	Policy # MI\ANTI\04\03\v03	Page 2 of 5
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual		

- 7. Place a cefoxitin disc in any available space remaining on the plate.
- 8. Incubate 35° C, in O₂ x 18-24 hours and record the zone diameters for the all cephalosporins as per NCCLS guidelines.

Figure 1. KB-ESBL Template



PROCEDURE MANUAL TORONTO MEDICAL LABORATORIES / MOUNT SINAI HOSPITAL MICROBIOLOGY DEPARTMENT

TML\MSH Microbiology Department	Policy # MI\ANTI\04\03\v03	Page 3 of 5
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual		

IV. Interpretation

Note: The following applies to cefpodoxime-nonsusceptible *E. coli*, *Klebsiella* species and *Proteus* species only.

- 1. Document zone size for all antibiotics.
- 2. Observe for **potentiation** of the inhibition zone (i.e. **increase** in the inhibition zone) of any one of cefpodoxime, ceftazidime, ceftriaxone or aztreonam when combined with clavulanic acid (enter **Y**es or **N**o to the "drug" named ESBL Inhibitor in the LIS).
- 3. If a reduction of zone of inhibition of any one of cefpodoxime, ceftazidime, ceftriaxone or aztreonam when combined with clavulanic acid is observed, recheck the identification of the isolate and repeat testing. Notify the charge technologist if result remains unchanged.

Class A ESBL present:

- i) Potentiation of the inhibition zone of any one of cefpodoxime, ceftazidime, ceftriaxone or aztreonam when combined with clavulanic acid (see below for examples of different patterns of potentiation that can be seen with organisms that contain Class A ESBLs)
- ii) Susceptibility to cefoxitin
- iii) Susceptibility or resistance to any one of ceftazidime, ceftriaxone or aztreonam



Class A and Class C ESBL present:

- i) Potentiation of the inhibition zone of any one of cefpodoxime, ceftazidime, ceftriaxone or aztreonam when combined with clavulanic acid
- ii) Resistant or Intermediate to cefoxitin.
- iii) Susceptibility or resistance to any one of ceftazidime, ceftriaxone or aztreonam

Class C-ESBL present:

- i) No potentiation with clavulanic acid
- ii) Resistance or Intermediate to cefoxitin
- iii) Resistance to any one of ceftazidime, ceftriaxone or aztreonam.

TML\MSH Microbiology Department	Policy # MI\ANTI\04\03\v03	Page 4 of 5
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual		

ESBL not Class A or Class C present:

- i) No potentiation with clavulanic acid
- ii) Susceptibility to cefoxitin
- iii) Resistance to any one of ceftazidime, ceftriaxone or aztreonam

ESBL absent:

- i) No potentiation with clavulanic acid
- ii) Susceptibility or resistance to cefoxitin
- iii) Susceptibilty to all of ceftazidime, ceftriaxone or aztreonam

V. Reporting

Reporting Comment	Potentiation of the inhibition zone of	Cefoxitin	Ceftazidime,
	any one of cefpodoxime, ceftazidime,		ceftriaxone
	ceftriaxone or aztreonam when		or aztreonam
	combined with clavulanic acid (enter		
	Y or N to the "drug" named ESBL		
	Inhibitor in the LIS)		
The susceptibility pattern suggests	Yes	S	S/R
that this organism contains a class			
A extended spectrum beta-			
lactamase (ESBL).			
The susceptibility pattern suggests	Yes	I/R	S/R
that this organism contains class A			
and C extended spectrum beta-			
lactamases (ESBL).			
The susceptibility pattern suggests	No	I/R	R
that this organism contains a class			
C extended spectrum beta-			
lactamase (ESBL).			
The susceptibility pattern suggests	No	S	R
that this organism contains an			
extended spectrum beta-lactamase			
(ESBL) other than class A or C.			
Not ESBL – no reporting comment	No	S/R	S

VII. <u>References</u>

1. Cormican MG, Marshall SA, Jones RN. 1996. Detection of extended-spectrum beta-lactamase (ESBL)-producing strains by the Etest ESBL screen. J. Clin. Microbiol. 8:1880-1884.

TML\MSH Microbiology Department	Policy # MI\ANTI\04\03\v03	Page 5 of 5
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual		

- 2. Jacoby GA, Han P. 1996. Detection of extended-spectrum beta-lactamases in clinical isolates of *Klebsiella pneumoniae* and *Escherichia coli*. J. Clin. Microbiol. 34:908-911.
- 3. Jacoby GA. 1994. Genetics of extended-spectrum beta-lactamases. Eur. J. Clin. Microbiol. Infect. Dis. Suppl 1,13:2-11.
- 4. Livermore DM. 1995. Beta-lactamases in laboratory and clinical resistance. Clin. Microbiol. Rev.8:557-584.
- 5. Moland ES, Sanders CC, Thomson KS. 1998. Can results obtained with commercially available MicroScan microdilution panels serve as an indictor of beta-lactamase (ESBL)-production among *Escherichia coli* and *Klebsiella* isolates with hidden resistance to expanded spectrum cephalosporins and aztreonam? J. Clin. Microbiol. 36:2575-2579.
- NCCLS. Performance standards for antimicrobial susceptibility testing; Eighth informational supplement. NCCLS document M100-SB [ISBN 1-26238-337-X]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 1998)
- Sanders CC, Barry AL, Washington JA, Shubert C, Moland ES, Traczewski MM, Knapp C, Mulder R. 1996. Detection of extended-spectrum beta-lactamase-producing members of the family Enterobacteriaceae with the Vitek ESBL test. J. Clin. Microbiol. 34:2997-3001.
- Swenson JM, Ferraro MJ, Jorgensen JH, Tenover FC. 1996. Can different breakpoints be used to detect extended-spectrum beta-lactamase (ESBL)-producing organisms? Abstr D-22, ICAAC meeting, New Orleans, Louisianna, Sept 15 -18, 1996.
- 9. Vercauteren E, Descheemaeker P, Ieven M, Sanders CC, Goossens H. 1997. Comparison of screening methods for detection of extended spectrum beta-lactamases and their prevalence among blood isolates of *Escherichia coli* and *Klebsiella* spp. in a Belgian teaching hospital. J. Clin. Microbiol. 35:2191-2197.

TML\MSH Microbiology Department	Policy # MI\ANTI\04\04\v03 Page	1 of 2
Policy & Procedure Manual		
Section: Antimicrobial Susceptibility Testing	Subject Title: Appendix IV – beta-lact	tamase
Manual	Testing	
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000	
Approved by: Laboratory Director	Revision Date: November 21, 2005	

APPENDIX IV - BETA - LACTAMASE TESTING

Introduction

This method can be used to detect beta-lactamase production in *Haemophilus* spp., *Neisseria gonorrhoeae*, *M. catarrhalis*, enterococci and staphylococci. Refer to <u>CRITERIA FOR</u> <u>SUSCEPTIBILITY TESTING</u> for when testing for beta-lactamase is required. The current method used in the laboratory is the chromogenic cephalosporin [cefinase (nitrocefin) disk] test. The test organism is inoculated directly onto the filter paper disk impregnated with nitrocefin. If the organism produces beta-lactamase, it will hydrolyze the chromogenic cephalosporin causing an electron shift that result in a coloured product.

II. <u>Materials</u>

Cefinase disks (BBL) (store refrigerated) Sterile distilled water Microscope slides Sterile Pasteur pipettes

III. Procedure

- 1. Using forceps remove the required number of disks from the dispenser and place on a microscope slide. Use 1 disk per organism.
- 2. Using a sterile Pasteur pipette, moisten each disk with a drop of sterile water.
- 3. With a sterile loop or applicator stick, pick several similar colonies from the agar plate and smear onto the surface of the disk.
- 4. Observe the disk for up to 5 minutes for a colour change. For staphylococci, observe the disks for up to 60 minutes.

IV. Interpretation

Positive:	the appearance of a pink colour
Negative:	no colour change

TML\MSH Microbiology Department	Policy # MI\ANTI\04\04\v03	Page 2 of 2
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual		

V. Quality Control

Set up positive and negative controls whenever a test is performed.

Positive:	<i>H. influenzae</i> ATCC 35056, β-lactamase positive.
Negative:	H. influenzae ATCC 10211, ß-lactamase negative.

VI. <u>Reference</u>

TML\MSH Microbiology Department Policy & Procedure Manual	Policy #MI\ANTI\04\05\v03	Page 1 of 2
Section: Antimicrobial Susceptibility Testing	Subject Title: Appendix V – Ox	acillin Screen
Manual	Plate	
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000	
Approved by: Laboratory Director	Revision Date: November 21, 20	005

APPENDIX V - OXACILLIN SCREEN PLATE

I. Introduction

This is an agar dilution method using a single concentration of oxacillin incorporated into Mueller Hinton (MH) agar to screen for resistant strains of *S. aureus*.

II. Materials

Control plate (MH with 4% NaCl) Screen plate (MH with 4% NaCl and 6 µg/mL OX) VITEK colorimeter Sterile saline Sterile swabs

III. Procedure

- 1. Using the VITEK colorimeter, prepare a suspension of the test organism (from solid medium after overnight culture) in sterile saline equivalent to a 0.5 McFarland standard using isolated colonies (inoculum prepared for VITEK can be used).
- 2. Using a sterile swab, spot inoculate the suspension onto the screen and control plates. Numerous organisms can be tested on one plate.
- 3. After the inocula have dried, incubate the plate at 35° C, O₂ for 24 hours.
- 4. All resistant isolates on the screen plate must be checked for purity (e.g. Gram stain, *S. aureus* tube coagulase or slide agglutination and sub-culture). The resistance must be confirmed by seting up a Denka MRSA Screen for confirmation of MRSA (Refer to DENKA MRSA SCREEN). Send a preliminary report as ISOLATE: Methicillin-resistant *S. aureus* and report to infection control.

IV. Interpretation

Growth on the screen plate indicates that the organism is methicillin resistant and therefore is considered resistant to <u>all</u> beta-lactam Antimicrobials (eg. penicillin, oxacillin, cephalosporins). Note: There must be growth on the control plate.

TML\MSH Microbiology Department	Policy #MI\ANTI\04\05\v03	Page 2 of 2
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual		

V. Quality Control

Controls must be tested each day. The organisms are to be sub-cultured from the TSB slant (in fridge) to Blood agar each day.

Sensitive:	S. aureus ATCC 29213
Resistant:	S. aureus LPTP 8610-1 S. aureus ATCC 43300

VI. Reference

Clinical and Laboratory Standards Institute (CLSI) Document – Methods for Dilution Antimicrobial Susceptibility Testing M7-A6, 2003.

TML/MSH Microbiology Department Policy & Procedure Manual	Policy #MI\ANTI\04\06\v02	Page 1 of 2
Section: Antimicrobial Susceptibility Testing	Subject Title: Appendix VI - De	enka MRSA
Manual	Screen	
Issued by: LABORATORY MANAGER	Original Date: July 31, 2000	
	\mathbf{D} '' \mathbf{D} () \mathbf{M} 1 01 0()0 <i>5</i>

APPENDIX VI - DENKA MRSA SCREEN

Principle

To be used as a screening test for the detection of Methicillin Resistant *S. aureus* (MRSA) from isolated colonies.

Reagents

MRSA screening kit Microcentrifuge tubes Boiling water bath Wooden sticks Loops Micropipettes

Method

- 1. Add 4 drops of extraction reagent #1 into a microcentrifuge tube.
- 2. Take one heavy loopful of the *Staphylococcus aureus* colonies from a blood plate and suspend the cells in the microcentrifuge tube.
- 3. Place in a boiling water bath for 3 minutes.
- 4. Remove microcentrifuge tube and let cool to room temperature.
- 5. Add one drop of extraction reagent #2 to the tube and mix well.
- 6. Centrifuge at high speed for 5 minutes.
- 7. For each specimen to be tested, allot and label one circle of the test card for testing with sensitized latex and one with control latex.
- 8. Place 50 microliters of the specimen onto 2 of the test circles and add one drop of the sensitized cells to one circle and one drop of the latex control to the other.
- 9. Mix the sample and latex together.
- 10. Rotate the card by hand for 3 minutes and observe for agglutination.

TML/MSH Microbiology Department	Policy #MI\ANTI\04\06\v02	Page 2 of 2
Policy & Procedure Manual		
Section: Antimicrobial Susceptibility Testing	Subject Title: Appendix VI - Denka	
Manual	MRSA Screen	

11. If negative rotate for another 3 minutes.

Interpretation

Sensitized latex	Control latex	Result
+	_	MRSA
-	-	Not MRSA
+ or -	+	Indeterminant

For DENKA – negative isolates:

- a. Perform Oxacillin-induced DENKA test:
 - i Sub-culture isolate onto a BA and add a Oxacillin disc. Incubate in O_2 at 35°C overnight.
 - ii Repeat DENKA using growth around the disc.
- b. If Oxacillin-induced DENKA is positive and it is a confirmed *S. aureus* report as MRSA.

If Oxacillin-induced DENKA is negative, check the Vitek Oxcillin MIC:

- a. Oxacillin =>4mcg/L, isolate is a BORSA; report as *S. aureus* with Isolate comment "This organism is resistant to cloxacillin by a mechanism different from that in typical MRSA. Consultation with a Microbiologist or Infectious Disease physician is advised."
- b. Oxacillin <4mcg/L, isolate is NOT a MRSA; report as *S. aureus* oxacillin-S.

Quality Control

Positive and negative controls must be set up once per week.

Positive:	S. aureus (ATCC 43330)
Negative:	S. aureus (ATCC 29213)

<u>Reference</u>

1. Denka Seiken Co., Ltd., Tokyo, Japan, Denka MRSA Screen package insert June 1998.

TML\MSH Microbiology Department Policy & Procedure Manual	Policy #MI\ANTI\04\07\v03	Page 1 of 3
Section: Antimicrobial Susceptibility Testing Manual	Subject Title: Appendix VII – S Bacteriostatic & Bactericidal T	erum Fitres
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000	
Approved by: Laboratory Director	Revision Date: November 21, 20	005

APPENDIX VII - SERUM BACTERIOSTATIC & BACTERICIDAL TITRES

I. <u>Introduction</u>

In the treatment of a patient with bacterial endocarditis or osteomyelitis, it is often important to know whether the prescribed dosages of Antimicrobials are achieving blood levels sufficiently high enough to kill the causative organism.

The bacteriostatic level is the dilution of serum that inhibits visible bacterial growth; the bactericidal level is the serum dilution that kills 99.9% of the initial inoculum.

NOTE: This test is to be performed only with the approval of a microbiologist.

II. <u>Specimen Collection</u>

The dose, the time the dose was given, and the time of collection must be recorded on the requisition. Pre- and post-dose blood specimens are obtained in serum separator tubes. The predose blood specimen is drawn immediately before administering the next dose of Antimicrobial in order to evaluate the pre (trough) level. Blood for the post-dose (peak) level should be drawn 1 hour after an intravenous infusion has been started, 1 hour after an intramuscular dose and 1 to 2 hours after an oral dose.

III. <u>Reagents/Materials/Media</u>

Mueller Hinton Broth (MHB) (100 mL) Blood Agar (BA) Sterile 13 x 100 mm glass tubes Sterile 1.0 mL pipettes Sterile yellow pipette tips Test tube racks Pipetter Precision pipette to deliver 20 ?L

IV. <u>Procedure</u>

A. Processing of Specimens

Upon arrival in the laboratory, centrifuge the blood and aseptically transfer the serum into a sterile vial.

TML\MSH Microbiology Department	Policy # MI\ANTI\04\07\v03	Page 2 of 3
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual		

B. Preparation of bacterial suspension:

Inoculate several colonies of a pure culture of the patient's organism (overnight subculture) into 5 mL MHB. Incubate on a shaker at 36° C for a minimum of 3 hours or until it achieves turbidity greater than the 0.5 McFarland standard (approximately 1.3 x 10^{8} CFU/mL).

- C. Serum dilution:
 - a. Place 12 sterile test tubes in a rack for each serum sample to be diluted.
 - b. Number the tubes 1 to 12.
 - c. Aseptically pipette 1.0 mL of patient's serum into tubes 1 and 2.
 - d. Aseptically pipette 1.0 mL of MHB into tubes 2 12.
 - e. With a new 1.0 mL sterile pipette transfer 1.0 mL of serum from tube 2 to tube 3. Mix well.
 - f. Serially dilute the serum by sequentially transferring 1.0 mL of the mixture through to tube 10. Discard 1.0 mL of the mixture from tube 10. No serum is to be added to tube 11 (positive inoculum control) or to tube 12 (broth sterility control). The final dilution of serum in tube 10 is 1:512 and final volume in all tubes should be 1.0 mL.
- D. Inoculating Broth
 - a. Using the Vitek colourimeter, dilute the bacterial suspension to the turbidity of the 0.5 McFarland standard using MHB.
 - b. Prepare a 1:4 dilution of the standardized inoculum by adding 1.0 mL of inoculum to 3.0 mL MHB. Mix well.
 - c. Using a precision pipette, dispense 20 μL (0.02 mL) of diluted inoculum into tubes 1 through 11. To inoculate, insert the pipette tip well under the surface of the Antimicrobial containing serum broth mixture. AVOID ANY CONTACT BETWEEN THE TIP AND THE WALLS OF THE TUBE to prevent transfer of organisms to the inside of tube above the meniscus. Mix by flushing 2 or 3 times without creating air bubbles or splashing. Use a new tip for each tube.
 - d. Incubate all tubes at 37° C for 20 hours in a CO₂-free incubator.
 - e. From the 1:4 dilution of the standardized inoculum, dilute 1:250 in MHB (0.1 mL in 24.9 mL MHB) to achieve an inoculum of 10^5 CFU/mL.
 - f. Perform a colony count to confirm the bacterial count in the final inoculum. Transfer 0.001 mL of diluted inoculum to BA by using a urine loop and distribute evenly on the surface of a BA plate.
 - g. Incubate the BA plate overnight at 35°C.

TML\MSH Microbiology Department	Policy # MI\ANTI\04\07\v03	Page 3 of 3
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual		

Determination of serum bacteriostatic titres

- 1. After incubation, tube 12 (broth sterility control) should be clear while tube 11 (positive inoculum control) should be turbid.
- 2. Record the colony count. The colony count plate should have 75-150 colonies. If the colony count is <75 or >150 consult the charge technologist before reading the tubes.
- 3. The highest dilution of serum that completely inhibits visible growth represents the bacteriostatic titre.

Determination of serum bactericidal titre

- 1. Vortex all tubes without visible growth for 15 seconds.
- 2. Use a urine loop to subculture all of the clear tubes onto 1/4 BA. Incubate at 37°C for 18 hours.
- 3. After incubation, read the plates and record the colony count.
- 4. The first dilution showing 99.9% killing activity (ie. no growth on sub-culture) is reported as the serum bactericidal titre.

V. <u>Reporting Results</u>

Telephone all results when available. Report as follows and give a copy of the report to the microbiologist:

Pre-dose serum bacteriostatic titre -Pre-dose serum bactericidal titre -

Post-dose serum bacteriostatic titre -Post-dose serum bactericidal titre -

VI. <u>Reference</u>

National Committee for Clinical Laboratory Standards. Methodology for the Serum Bactericidal Test, NCCLS Document M21-P,Vol. 7, No. 1, 1987.

TML\MSH Microbiology Department	Policy # MI\ANTI\04\08\v02	Page 1 of 1
Policy & Procedure Manual		
Section: Antimicrobial Susceptibility Testing	Subject Title: Appendix VIII –	Broth
Manual	Macrodilution and Agar	
	Dilution	
Issued by: LABORATORY MANAGER	Dilution Original Date: January 10, 2000	
Issued by: LABORATORY MANAGER Approved by: Laboratory Director	Dilution Original Date: January 10, 2000 Revision Date: November 21, 20	005

APPENDIX VIII - BROTH MACRODILUTION AND AGAR DILUTION

I. Introduction

These tests are not routinely done and will only be performed following consultation with a microbiologist. Refer to the NCCLS standard M7-A6 (January 2003) for methodology.

TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI\ANTI\04\09\v03	Page 1 of 2
Section: Antimicrobial Susceptibility Testing	Subject Title: Appendix IX – Br	oth
Manual	Microdilution MI	C
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000	
Approved by: Laboratory Director	Revision Date: November 21, 20)05

APPEDIX IX - BROTH MICRODILUTION MIC

I. Introduction

Dilution susceptibility testing methods are used to determine the minimal concentration of an antimicrobial agent required to inhibit or kill a microorganism. Antimicrobial agents are usually tested at log₂ (twofold) serial dilutions, and the lowest concentration that inhibits visible growth of an organism is regarded as the MIC. The concentration range used may vary with the drug, the organism tested, and the site of infection. The method and principles of the microdilution method is essentially the same as the macrodilution method except that the antimicrobial dilutions are in 0.1 ml volumes contained in wells of a microdilution tray (usually 96 well trays). Results obtained may be reported as the actual MIC or categorically as Susceptible (S), Intermediate (I), or Resistant (R). Interpretive categories are published and up-dated regularly by NCCLS.

II. <u>Materials</u>

Sterile saline Transfer pipettes Sterile distilled water Vitek colorimeter MIC microtitre panel Inoculator (tray and lid)

III. Procedure

- 1. Prepare a suspension of the test organism in sterile saline equivalent to a 0.5 McFarland standard using isolated colonies.
- 2. Transfer 1.5 mL of the suspension to the bottom of the inoculating tray and add approximately 40 mL of sterile distilled water.
- 3. Aseptically replace the transfer lid into the inoculating tray making sure no bubbles are under the prongs.
- 4. Lift the transfer lid and center it over the previously thawed MIC panel.

TML\MSH Microbiology Department	Policy # MI\ANTI\04\08\v03	Page 2 of 2
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual		

- 5. Lower the transfer lid into the panel so the prongs enter all wells except the control well.
- 6. Remove transfer lid and cover the panel with a dummy MIC panel.
- 7. Incubate for 18-24 hours at 35°C, depending on drugs/bugs tested.

IV. <u>Interpretation</u>

The highest dilution of the Antimicrobial that completely inhibits visible growth represents the minimum inhibitory concentration (MIC).

V. Quality Control

Panels are Quality Controlled with the appropriate ATCC control organisms.

VI. <u>Reference</u>

Clinical and Laboratory Standards Institute (CLSI) Document – Methods for Dilution Antimicrobial Susceptibility Testing M7-A6, 2003.

TML\MSH Microbiology Department	Policy ## MI\ANTI\04\10\v02	Page 1 of 2
Policy & Procedure Manual		
Section: Antimicrobial Susceptibility Testing	Subject Title: Appendix X – e-te	est
Manual		
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000	
Approved by: Laboratory Director	Revision Date: November 21, 20)05

APPENDIX X - E-TEST

I. Introduction

The E test (also known as the Gradinet Diffusion Method) is based on the same principle as the disk diffusion method. It is an in vitro method for quantitative antimicrobial susceptibility testing whereby a preformed antimicrobial gradient from a plastic-coated strip diffuses into an agar medium inoculated with the test organism. The MIC is read directly from a scale on the top of the strip at a point where the ellipse of organism growth inhibition intercepts the strip.

II. <u>Materials</u>

E-test strips (store frozen) Mueller Hinton Agar (MH) Mueller Hinton Blood Agar (MHBA) Haemophilus Test Media (HTM) Trypticase Soy Broth (TSB) (3 mL) VITEK colorimeter Sterile saline Sterile swabs

III. Procedures

- 1. Allow E-test strips to come to room temperature before opening the container.
- Using the Vitek colorimeter, prepare a suspension of the test organism in sterile saline equivalent to a 0.5 McFarland standard using isolated colonies. If there is not enough growth, inoculate the organism into TSB, and incubate at 35°C for 2-4 hours or until it reaches the turbidity of a 0.5 McFarland standard.
 For mucoid organisms, adjust suspension to 1 McFarland standard.
- 3. Using a sterile cotton swab, inoculate the organism onto an appropriate agar plate, streaking in 3 directions over the entire agar surface. For organisms that grow rapidly use MH agar. For *Haemophilus* species use HTM and for *S. pneumoniae* and viridans streptococci use MHBA. For other organisms that do not grow on MH, use MHBA.
- 4. Using forceps, apply the appropriate Antimicrobial strips onto the agar. Use one plate per Antimicrobial strip.

TML\MSH Microbiology Department	Policy # MI\ANTI\04\10\v02	Page 2 of 2
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual		

5. Incubate plates as follows:

Haemophilus species - CO_2 , 35°C x 18 hours *Streptococcus pneumoniae* - CO_2 , 35°C x 20 to 24 hours *Staphylococcus aureus* and *Enterococcus* species for Methicillin and Vancomycin - O_2 , 35°C x 24 hours Others - O_2 , 35°C x 18 hours

IV. Interpretation

After incubation read the MIC value at the point of intersection between the zone edge and the E-test strip. Since E-test comprises a continuous gradient, MIC values in between two-fold dilutions can be obtained. Always round up these values to the next two-fold dilution before interpretation.

See e-test Reading Guide

For polymyxin B see <u>etest Polymyxin B Reading</u>

V. Quality Control

The 3 E-test strips are tested once weekly with *S. aureus* ATCC 29213. The organism is subcultured from the TSB slant (in fridge) to BA the day before setting up the QC.

Expected Results^{*}:

MIC

0.25-2.0 μg/mL
4.0-16.0 μg/mL
1.0-8.0 µg/mL
1.0-4.0 ?g/mL

? As per NCCLS document M100-S15 (M7), Table 3, January 2005.

VI. <u>Reference</u>

TML\MSH Microbiology Department	Policy # MI\ANTI\04\11\v04 Page 1 of 2	
Policy & Procedure Manual		
Section: Antimicrobial Susceptibility Testing	Subject Title: Appendix XI – Vancomycin	
Manual	and High Level Aminoglycoside	
	Testing	
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000	
Approved by: Laboratory Director	Revision Date: November 21, 2005	

APPENDIX XI - VANCOMYCIN & HIGH LEVEL AMINOGLYCOSIDE TESTING

I. Introduction

Synergy between ampicillin, pencillin or vancomycin and an aminoglycoside for *Enterococcus* species can be predicted by high level aminoglycoside (HLAR - gentamicin and streptomycin) screening test. Vancomycin resistance of *Enterococcus* species can be detected by BHI vancomycin agar screen plate containing 6mg/L of vancomycin.

II. <u>Materials</u>

Control plate (Brain Heart Infusion Agar) Entero HLAR and Vancomycin Screen plates VITEK colorimeter Sterile saline Sterile swabs

III. Procedure

- 1. Using the VITEK colorimeter, prepare a 0.5 McFarland suspension in sterile saline (inoculum from VITEK can be used).
- 2. Using a sterile swab, spot inoculate the suspension onto each of the test and control plates.
- 3. After the inocula have dried, incubate the plate at 35°C for up to 48 hours.

IV. Interpretation

Check the control plate for adequate growth. Then check the drug plates for absence or presence of growth; any growth is considered significant. Read plates at 24 hours and record results. If there is no growth on the drug plates, re-incubate plate for an additional 24 hours.

TML\MSH Microbiology Department	Policy # MI\ANTI\04\11\v04	Page 2 of 2
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual		

V. Quality Control

Control strains are tested with each plate.

	Expected results of each quadrant				
	C	G	S	V	
E. faecalis (ATCC 19966)	+	+	-	-	
E. faecalis (ATCC 49533)	+	-	+	-	
E. gallinarum (ATCC 35038)	+	-	-	+	
E. casseliflavus (ATCC 12755)	+	-	+	+	
S. aureus (ATCC 29213)	+			-	

C = Growth Control; G = Gentamicin; S = Streptomycin; V = Vancomycin

VI. <u>Reporting Results</u>

Blood cultures and sterile sites:

- ? If high level gentamicin is **susceptible** (regardless of streptomycin result) report as:
 - Serious enterococcal infections may require an aminoglycoside for synergy.
 Please contact the Medical Microbiologist for treatment advice".
- ? If high level gentamicin is resistant (regardless of streptomycin result) report as:

 « "This organism is high level aminoglycoside resistant. Please contact the Medical Microbiologist for treatment advice".
- ? Record the streptomycin result in the LIS. Report result only upon request.

Urines and other sites:

- ? Do not report HLAR.
- ? Report Vancomycin result as per specimen type specific reporting tables.

VII. <u>Reference</u>

PML Technical Manual data sheet No. 323, Nov. 1993.

TML\MSH Microbiology Department Policy & Procedure Manual	Policy #MI\ANTI\04\12\v02	Page 1 of 2
Section: Antimicrobial Susceptibility Testing	Subject Title: Appendix XII – V	ancomycin
Manual	Screen Plate for St	taphylococcus
		1 v
Issued by: LABORATORY MANAGER	Original Date: May 31, 2004	
Issued by: LABORATORY MANAGER Approved by: Laboratory Director	Original Date: May 31, 2004 Revision Date: November 21, 20	005

APPENDIX XII - VANCOMYCIN SCREEN PLATE FOR STAPHYLOCOCCUS

I. Introduction

This is an agar dilution method using a single concentration of Vancomycin ($6 \mu g/mL$) incorporated into Brain Heart Infusion (BHI) agar to screen for resistant strains of *Staphylococcus aureus*.

II. <u>Materials</u>

Control plate (Brain Heart Infusion Agar) Screen plate (Brain Heart Infusion Agar with 6 µg/mL Vancomycin) VITEK colorimeter Sterile saline Sterile swabs

III. Procedure

- 1. Using the VITEK colorimeter, prepare a 0.5 McFarland suspension in sterile saline (inoculum from VITEK can be used).
- 2. Using a sterile swab, spot inoculate the suspension onto each of the test and control plates.
- 3. After the inocula have dried, incubate the plate at 35° C in O₂ for 24 hours.
- 4. All resistant isolates on the screen plate must be checked for purity (e.g. Gram stain, tube coagulase or slide agglutination and sub-culture). The resistance must be confirmed by repeating the screen plate and setting up e-test (Refer to APPENDIX X E-TEST). Report preliminary results to infection control.

IV. Interpretation

Check the control plate for adequate growth. Then check the drug plates for absence or presence of growth; any growth is considered significant. Read plates at 24 hours and record results. Note: There must be growth on the control plate.

TML\MSH Microbiology Department Policy & Procedure Manual	Policy #MI\ANTI\04\12\v02	Page 2 of 2
Antimicrobial Susceptibility Testing Manual		

V. Quality Control

Control strains are tested on every plate.

	С	V
E. faecalis (ATCC 19966)	+	_
E. faecalis (ATCC 49533)	+	-
E. gallinarum (ATCC 35038)	+	+
E. casseliflavus (ATCC 12755)	+	+
S. aureus (ATCC 29213)	+	-

C = Growth Control; V = Vancomycin

VI. <u>Reference</u>

bioMerieux procedure update May 28, 2004

TML\MSH Microbiology Department	Policy #MI\ANTI\04\13\v03	Page 1 of 6
Policy & Procedure Manual		
Section: Antimicrobial Susceptibility Testing	Subject Title: Appendix XIII -	_
Manual	Antimicrobial A	bbreviations
Issued by: LABORATORY MANAGER	Original Date: January 10, 200	00
Issued by: LABORATORY MANAGER Approved by: Laboratory Director	Original Date: January 10, 200 Revision Date: November 21,	00 2005

APPENDIX XIII – Antimicrobial Abbreviations

Abbrevaitions - Antimicrobial Disks

ANTIMICROBIAL	DISK (Manufacturer)	Concentration (?g)
Amikacin	AK (Oxoid)	30
Amoxacillin/Clavulanic Acid	AMC	30
Ampicillin	AMP (Oxoid)	10
Aztreonam	ATM	30
Cefazolin	KZ (Oxoid)	30
Cefepime	FEP	30
Cefixime	CFM	5
Cefotaxime	СТХ	30
Cefotetan	CTT (Gen. Diag.)	30
Cefoxitin	FOX (Oxoid)	30
Ceftazidime	CAZ (Oxoid)	30
Ceftriaxone	CRO (Oxoid)	30
Cefuroxime	CXM	30
Cephalothin	KF	30
Cefpodoxime	CPD	10
Chloramphenicol	С	30
Ciprofloxacin	CIP (Oxoid)	5
Clarithromycin	CLR	15
Clindamycin	DA (Oxoid)	2
Colistin	СТ	10
Cotrimoxazole	SXT (Oxoid)	
Erythromycin	E (Oxoid)	15
Fusidic Acid	FD	10
Gentamicin	CN (Oxoid)	10
Imipenem	IPM (Difco)	10
Levofloxacin	LVX	5
Meropenem	MEM	10
Metronidazole	MTZ (Oxoid)	
Minocycline	MH	30
Mupirocin	MUP	5
Nalidixic Acid	NA	30
Nitrofurantoin	F (Oxoid)	300
Norfloxacin	NOR (BBL or Difco)	10
Novobiocin	NV	5

PROCEDURE MANUAL

TML\MSH Microbiology Department	Policy #MI\ANTI\04\13\v03	Page 2 of 6
Policy & Procedure Manual		
Antimicrobial Sus ceptibility Testing Manual		

ANTIMICROBIAL	DISK (Manufacturer)	Concentration (?g)
Oxacillin	OX (Oxoid)	1
Penicillin	P (Oxoid)	10
Piperacillin	PRL (Oxoid)	100
Pipercillin/Tazobactam	TZP	110
Rifampin	RA	5
Teicoplanin	TE	30
Tetracycline	TE (Oxoid)	30
Timentin	TIM	85
Tobramycin	TOB (Oxoid)	10
Trimethoprim/Sulfamethoxazole	SXT	
Vancomycin	VA (Oxoid)	30

TML\MSH Microbiology Department	Policy #MI\ANTI\04\13\v03	Page 3 of 6
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual		

Abbreviations – E-Test Strips

ANTIMICROBIAL	ABBREVIATION	CONCENTRATION (?g)
Amikacin	AK	30
Ampicillin	AM	30
Azithromycin	AZ	30
Cefotaxime	CT	100
Cefotetan	CN	100
Cefoxitin	FX	100
Ceftazidime	TZ	30
Ceftriaxone (For CHC Pneumo)	TX	30
Cefuroxime	XM	30
Cephalotin	CE	30
Chloremphenicol	CL	30
Ciprofloxacin	CI	30
Clarithromycin	СН	30
Clindamycin	СМ	100
Colistin	CO	
Doxycycline	DC	30
Erythromycin	EM	30
Fusidic Acid		30
Gentamycin (Low Level)	GM	30
High Level Gentamicin	GM	30
Imipenem	IP	30
Metronidazole	MZ	100
Mupirocin	MU	100
Oxacillin	OX	30
Penicillin	PG	100
Pipercillin	PP	100
Pipercillin/Tazobactam	PTC	30
Polymyxin B	PO	
Streptomycin (High Level)	SM	30
Teicoplanin		30
Tetracycline		30
Vancomycin	VA	100

TML\MSH Microbiology Department Policy & Procedure Manual Policy #MI\ANTI\04\13\v03

Page 4 of 6

Antimicrobial Susceptibility Testing Manual

Abbreviations – LIS (Soft Computer Corporation)

ANTIMICROBIAL	ABBREVIATION
Amikacin	an
Amoxicillin	amx
Amoxicillin / Clavalanic Acid	amc
Ampicillin	am
Ampicillin / Sulbactam	ama
Azithromicin	azi
Aztreonam	azm
?eta-lactamase	blac
?-lactamase	beta
Carbenicillin	cb
Cefaclor	ccl
Cefamandole	cm
Cefazolin	CZ
Cefipime	сро
Cefixime	cfm
Cefotaxime	tax
Cefotetan	cte
Cefoxitin	fox
Cefpodoxime	cpd
Cefpodoxime	cpod
Cefpodoxime / Clavulanic Acid	cpodc
Ceftazidime	taz
Ceftizoxime	ZOX
Ceftriaxone	ctr
Cefuroxime	roxh
Cefuroxime – Axetil	roxa
Cefuroxime-sodium	rox
Cephalothin	cf
Chloramphenicol	С
Ciprofloxacin	cip
Clarifloxacin	clar
Clarithromycin	cla
Clinafloxacin	cflox
Clindamycin	сс
Cloxacillin	clx
Colistin	ct
Dalfopristin	dalfo
Doxycycline	dx
D-zone	dzone

TML\MSH Microbiology Department	Policy #MI\ANTI\04\13\v03	Page 5 of 6
Policy & Procedure Manual		
Antimicrobial Susceptibility Testing Manual		

ANTIMICROBIAL	ABBREVIATION
Erythromycin	e
ESBL Potentiation	esbinh
Everninomycin	ever
Fusidic Acid	fa
Gatifloxacin	gat
Gentamicin	gm
Gentamicin 2000	gm2000
Gentamicin 500	gm500
Imipenem	imi
Kanamycin	k
Levofloxacin	lev
Linezonomycin	linezo
Meropenem	mem
Metromidazole	mtz
Mezlocillin	mz
Minocycline	mn
Mupirocin	mup
Nalidixic Acid	na
Netilmicin	net
Nitrofurantoin	fd
Norfloxacin	nor
Ofloxacin	ofx
Oxacillin	OX
Penicillin	peng
Piperacillin	pip
Pipercillin / Tazobactam	tzp
Polymyxin B	pb
Pristinamycin	pris
Ramoplanin	ramo
Rifampin	rif
Streptomycin	strep
Streptomycin 2000	st2000
Sulfisoxazole	SOX
Synercid	suncd
Teicoplanin	tei
Tetracycline	tet
Ticarcillin	tic
Ticarcillin/Clavulanic Acid	tcc

TML\MSH Microbiology Department Policy & Procedure Manual	Policy #MI\ANTI\04\13\v03	Page 6 of 6
Antimicrobial Susceptibility Testing Manual		

ANTIMICROBIAL	ABBREVIATION
Tobramycin	tob
Trimethoprim	tmp
Trimothoprim/sulfamethoxazole	sxt
Vancomycin	va

OXACILLIN AND VANCOMYCIN SCREEN RECORDING SHEET FOR S. aureus Policy # MI\ANTI\04\03a\v03

D	ate:	Set-up	by:		Lot	# (OXA	(VA	NC)_		Rea	d by:	18h 24h_		-	00	
#	ATCC Control Or Lab No.	OXA 18 h 24	VA 181	NC 1 24	#	Lab No.	O2 18 1	XA 1 24	VA 181	NC 1 24	#	Lab No.	OXA 2	18 h 4	VA 18 h	NC 1 24
1	OXA: 8610 VANC: 49532	R	s		14						27					
2	OXA: 43300 VANC: 49573	R	R		15						28					
3	OXA: 43387 VANC: 49533	S	S		16						29					
4	OXA: 29213 VANC: 29213	S	S		17						30					
5					18						31					
6					19						32					
7					20						33					
8					21						34					
9					22						35					
10					23						36					
11					24						37					
12					25											
13					26											

QUAD SCREEN RECORDING SHEET FOR ENTEROCOCCUS and CNST

	Data			Sot ur	- h			Lot#(vene):						Policy# MI	ANT	1\04	\08a\v	03	
	48h			Set-ul	5 Uy.			Lot# (valie)			(nL)	1)		Keau 0y. 1011		. 24	II	_	
#	Lab No.	V 181	'A h 24	GM 24h	S 24	Г R h 48	#	Lab No.	V 18 1	'A h 24	GM 24h	STR 24 h 48	#	Lab No.	VA 18 h 24		GM 24h	ST 24 h	Г R h 48
1	<i>E. faecalis</i> ATCC 49532	S		R	S		1 4						2 7						
2	E. gallinarum ATCC 49573	R		S	S		1 5						2 8						
3	E. faecalis ATCC 49533	s		S	R		1 6						2 9						
4	S .aureus ATCC 29213	S			-		1 7						3 0						
5							1 8						3 1						
6							1 9						3 2						
7							2 0						3 3						
8							2 1						3 4						
9							2 2						3 5						
1 0							2 3						3 6						
1							2 4						3 7						
1 2							2 5						3 8						
1 3							2 6						3 9						

TML\MSH Microbiology Department Policy & Procedure Manual	Policy #MI\ANTI\04\14\v03	Page 1 of 3				
Section: Antibiotic Susceptibility Testing Manual	Subject Title: Appendix XIV - VITEK					
Issued by: LABORATORY MANAGER	Original Date: January 10, 2000					
Approved by: Laboratory Director	Revision Date: February 5, 2002	· · ·				

APPENDIX XIV - VITEK

I. Introduction

The VITEK System is an automated, short-incubation micro broth dilution system capable of performing susceptibility testing of most rapidly growing gram-positive and gram-negative aerobic bacteria. It can yield a result in a period of 4 to 10 hours. Results are automatically transferred from the VITEK System to the Laboratory Information System (LIS) via a computer interface. Results are reported qualitatively as either Susceptible (S), Intermediate (I), or Resistant (R).

II. <u>Materials</u>

Cards (stored refrigerated):

GNI --- Gram negative identification card GPI --- Gram positive identification card GNS – 606 --- Gram negative susceptibility card GPS –105 --- Gram positive susceptibility card

1.8 mL 0.45% sterile saline VITEK colorimeter MLA pipette Pipette tips Transfer tubes Filling Stand

III. Procedure

A. Preparation of Inoculum:

Use fresh 24 hours culture.

If insufficient growth for setting up inoculum, sub-culture isolate onto a Blood Agar (BA) plate, incubate for a few hours and set up from a sweep of the growth.

Sub-culture frozen or freeze dried isolates twice before setting up any Vitek card.

TML\MSH Microbiology Department Policy & Procedure Manual	Policy # MI\ANTI\04\14\v03	Page 2 of 3
Antibiotic Susceptibility Testing Manual		

Gram Negative(GNI / GNS)Organism ?1.8 mL saline (1 McFarland) - 50 ?L ?1.8 mL salineGNIGNS

Gram Positive(GPI / GPS)Organism ?1.8 mL saline (0.5 McFarland) - 200 ? L ?1.8 mL salineGPIGPS

Vortex suspensions for ID card inoculations Inoculate cards within 15 minutes.

- B. Card Inoculation:
 - 1. Allow card to reach room temperature.
 - 2. Remove protective package liner, visually inspect for holes or cracks in the foil material.
 - 3. Label card with Last 6 digits of Lab # + isolate #. Isolate # should be identical in corresponding ID and susceptibility cards.
 - Mark external tests marks (positives) on appropriate cards: oxidase/catalase beta hemolysis/coagulase Corresponding ID and susceptibility cards should have identical external tests marked.
 - 5. Insert transfer tube into inlet port.
 - 6. Place inoculated test tube in filling stand.
 - 7. Place assembled card and transfer tube on filling stand with the long portion of the transfer tube into the test tube.
 - 8. Place filling stand with the assembled test kit into filling module.
 - 9. Press 'FILL' button once. 'READY' light will go out after a few moments.
 - 10. When the 'READY' light returns, remove filling stand from filling module.
 - 11. Cut transfer tube and seal card in heater/sealer module.
 - 12. Check card for bubbles.
 - 13. For ID Cards, remove bubbles by shaking.
 - 14. For sensitivity cards, DO NOT SHAKE. Small bubbles will not cause problems. If very large bubbles are present, discard the card and repeat.
 - 15. Insert card into tray.

16. Put filled card into reader/incubator within 15 minutes.

- 15. Check time and location of the next tray to be read by:
 - a) click on "**VITEK**" in main menu
 - b) click on "**READER**"
 - c) click on '**STATUS**" Open reader/incubator.

TML\MSH Microbiology Department	Policy # MI\ANTI\04\14\v03	Page 3 of 3
Policy & Procedure Manual		
Antibiotic Susceptibility Testing Manual		

Put new card into tray following 'next tray to be read'. Example: Next tray to be read: tray a put new card into tray b. (If b is full, trays c or d can also be used).

To Exit "STATUS" WINDOW:

- a) click on '**FILE**"
- b) click on 'QUIT"

NHI, ANI CARDS:

After inoculation, incubate in non-CO₂ incubator for 4 hours.

At amscommand, type **'nhi'** (not capitals) Enter ID # Catalase MTM Indole Results of card per NIH insert

C. Quality Control:

QC all new lots of **identification cards** once when received in the lab. No further QC is required if all results are within limits. Perform QCs for 5 consecutive days if out-of-control results are observed (see QA technologist).

QC all new lots of **susceptibility cards** once when received and once weekly when lot is in use. Perform daily QCs when any out-of-control results are observed (see QA technologist).

See VITEK manual for QC organisms to be set up and set up procedure.