

TML/MSH Microbiology Department Policy & Procedure Manual	Policy # MI/EDU/04/04/v01	Page 1 of 3
Section: Education Manual	Subject Title: Appendix IV - Tests	
Issued by: LABORATORY MANAGER	Original Date: June 18, 2001	
Approved by: Laboratory Director	Revision Date:	

APPENDIX IV - TESTS

1. Catalase
 - bacteria that contain cytochrome enzymes are catalase positive and those that don't are catalase negative
 - *Staphylococcus* +
 - *Streptococcus* and *Enterococcus* -
 - flood bacteria with 3% hydrogen peroxide and observe for bubbles catalase
 - a) H_2O_2 (3%) \rightarrow Catalase peroxide + H_2O
 - b) Catalase peroxide + $\text{H}_2\text{O}_2 \rightarrow$ Catalase $\text{H}_2\text{O} + \text{O}_2$

2. Coagulase
 - the ability to clot plasma
 - two different coagulase tests can be performed, a tube test for free coagulase and a slide test for bound coagulase, or clumping factor
 - *Staphylococcus aureus* +
 - Coagulase negative *Staphylococcus* -

3. Bile esculin
 - to detect beta glucoside which breaks down esculin to form a black precipitate due to the presence of ferric ions
 - *Enterococcus faecalis* +
 - Beta-haemolytic group B streptococcus -

Esculin (β glucoside) \rightarrow Esculetin + ferric ions (ferric citrate in medium)

↓

black precipitate

4. MUG
 - this test is mainly done on the urine bench
 - to detect beta-glucuronidase
 - colourless substrate is broken down to produce a yellow compound
 - Kovac's reagent is added to detect indole production
 - *Escherichia coli* is MUG + and indole +

TML/MSH Microbiology Department Policy & Procedure Manual	Policy # MI/EDU/04/04/v01	Page 2 of 3
Education Manual		

tryptophanase (produced by bacteria)
Tryptophan → ammonia + pyruvic acid + indole
↓
reacts with p-dimethylamino benzaldehyde (Kovac's reagent)
↓
quinoidal red violet

[Indole spot test uses Erlich's reagent 1% dimethylaminocinnamaldehyde]

↓
blue colour

5. Oxidase
- to test for the production of oxidase
 - spot inoculate organism on to a filter paper soaked with 1% tetramethylphenylene diamine dihydrochloride - positive is purple, negative is yellow
 - *Pseudomonas aeruginosa* +
 - *Escherichia coli* -

Oxidizing reaction

Reagent 1% Dimethyl or Tetramethyl para-phenylenediamine

↓ on colonies

↓

Indophenoloxidase (produced by bacteria)

↓

↓

Indophenol = black colonies (with dimethyl)

= magenta colonies (with tetramethyl)

The following tests are done on the stool bench to screen for *Salmonella*, *Shigella* and *Yersinia*.

6. Urease
- Detects urease production
 - Peptones in media utilized producing an alkalinity. Phenol red indicator
- Urea →→ ammonia

TML/MSH Microbiology Department Policy & Procedure Manual	Policy # MI/EDU/04/04/v01	Page 3 of 3
Education Manual		

7. TSI
- Triple sugar iron agar. 0.1% glucose, 1% sucrose and lactose.
 - If the glucose is fermented, only a small amount of acid will be produced, which will be neutralized by alkali from peptone metabolism along surface of slant.
 - Oxidation of peptone cannot take place in the anaerobic conditions in the depth of the medium.
 - Therefore when glucose only is fermented, the butt of the medium becomes yellow and the slant remains red.
 - If lactose or sucrose is fermented, the amount of acid produced is large enough to offset alkali production and a yellow slant is produced.
 - Production of hydrogen sulphide is shown by formation of iron sulphide from the ferrous sulphate.
8. ONPG
- To detect enzyme β -D-galactosidase in lactose fermenting organisms. O-nitrophenol-B-D-galactopyranoside $\rightarrow\rightarrow$ O-nitrophenol (yellow)
9. PPA
- To detect Phenylalanine deaminase production. Phenylalanine $\rightarrow\rightarrow$ Phenylpyruvic acid + FeCl_3 (ferric ions) = blue green