

MSH/TML Shared Microbiology Service Policy & Procedure Manual	Policy # M\PAR\05\13\v01	Page 1 of 2
Section: Parasitology Manual	Subject Title: Laboratory Procedures for Stool Examination	
Issued by: LABORATORY MANAGER	Original Date: March 13, 2000	
Approved by: Laboratory Director	Revision Date:	

Charcoal Culture for *Strongyloides*

PRINCIPLE

Culture techniques are useful to detect light numbers of larvae not found during normal concentration procedures. They also allow for easier speciation between strongyloides and hookworm due to development of rhabidiform larvae of hookworm. Finally culture techniques allow larvae to develop to the filariform stage to further aid in diagnosis. This method is less sensitive than the agar plate method since it relies on directly observing the larvae. Charcoal is used to maintain pH and to provide a medium in which the larvae can develop.

SPECIMEN

Fresh stool sample that has not been refrigerated.

SAFETY

- Assume that the sample contains filariform larvae and wear gloves and take measures to prevent the larvae from migrating out of the dish.
- Warn co-workers about the nature of your work.

PROCEDURE

This is considered to be a non-routine procedure therefore it should only be performed by experienced personnel.

1. Mix 20-40 g of fecal material in tap water until a thick suspension is obtained.
2. Add this mixture to a suitable covered dish to which has been added a similar volume of #10 granulated hardwood charcoal.
3. Mix the suspension well using a tongue depressor. Add enough water so that the charcoal will glisten but not enough to allow it to lay at the bottom of the dish. Keep at room temperature.

MSH/TML Shared Microbiology Service Policy & Procedure Manual	Policy # MI\PAR\05\13\v01	Page 2 of 2
Parasitology Manual		

4. Check culture daily for sufficient moisture and sprinkle water on the surface, if needed. *As infective filariform larvae could be present at any time during this procedure, caution must be used when any handling of the culture occurs.*
5. After 5-6 days the larvae, if present, can be harvested by using the Baermann technique.