

TML/MSH Microbiology Department Policy & Procedure Manual	<b>Policy # MI/VIR/16/05/v02</b>	Page 1 of 4
Section: <b>Virology Manual</b>	Subject Title: <b>Appendix V Direct Immunofluorescent Antibody (DFA) staining for Viral Culture Confirmation</b>	
Issued by: <b>LABORATORY MANAGER</b>	Original Date: March 14, 2001	
Approved by: Laboratory Director	Revision Date: October 10, 2003	

## Appendix V

### DIRECT IMMUNOFLUORESCENT ANTIBODY (DFA) STAINING

#### I. Introduction

The DFA staining technique is used to detect viruses either directly in patient specimens or which have been isolated in shell vial or tube cultures. The method consists of a single staining step using a virus-specific antibody which is conjugated with a fluorochrome. Viruses which we currently identify by DFA staining include HSV-1, HSV-2, VZV, CMV (late antigen) and respiratory viruses (SimulFluor stains for respiratory syncytial virus, parainfluenza, influenza, adenovirus).

#### II. Reagents and Materials

FITC-conjugated virus-specific antibody  
FITC/Rodamine-conjugated virus-specific antibody (SimulFluor)  
Phosphate Buffered Saline (PBS)  
dH<sub>2</sub>O  
cold acetone (4°C)  
mounting fluid  
sterile pipettes  
cytospin and accessories (for tube culture)  
humidified chamber  
glass slides  
coverslips  
paper towels for blotting

#### III. Procedure

##### 1. **Shell Vial**

This procedure is for staining of cells directly in shell vial. If staining a cytospin slide or slide made directly from a patient specimen, follow the tube culture procedure below.

TML/MSH Microbiology Department Policy & Procedure Manual	<b>Policy # MI/VIR/16/05/v02</b>	Page 2 of 4
<b>Virology Manual</b>		

- i. Discard cap. Remove maintenance medium from the shell vial using sterile pipette.
- ii. Add 1 mL of cold acetone. Cover with tray lid and let sit for 10 minutes.
- iii. Decant acetone and blot shell vial on paper towel.
- iv. Gently rinse with PBS from squirt bottle, filling vial 3/4 full. Decant PBS.
- v. Add 75µl (2 drops from bottle) of appropriate FITC-conjugated virus-specific antibody. Cover with tray lid.
- vi. Incubate at 36°C for 30 minutes.
- vii. Gently rinse with PBS from squirt bottle, filling vial 3/4 full. Decant PBS. Repeat.
- viii. Remove the coverslip from each shell vial and place cell side down onto a drop of mounting fluid on a glass slide.
- ix. For HSV 1, HSV 2, VZ and CMV, read using fluorescence microscope with the FITC/Evans Blue filter and the 40x objective.
- x. For respiratory viruses, read using fluorescence microscope with the FITC/Evans Blue filter and the 40x objective.

## **2. Tube Culture**

- i. Prepare cytospin slide from cell culture tube as outlined in Appendix XX.
- ii. Fix slide in cold acetone for 10 minutes in a coplin jar. Remove slide and air dry.
- iii. Add 20µl of appropriate FITC-conjugated antibody onto the fixed cytospin slide.
- iii. Incubate in a humidified chamber at 36°C for 30 minutes.
- iv. Wash each slide 3 times with fresh PBS for 2 minutes each in a coplin jar.
- v. Wash with distilled water for 1 minute in a coplin jar.
- vi. Wipe excess water from the slide without touching the cytospin preparation.

TML/MSH Microbiology Department Policy & Procedure Manual	<b>Policy # MI/VIR/16/05/v02</b>	Page 3 of 4
<b>Virology Manual</b>		

- vii. Mount using coverslip and mounting fluid.
- xi. For HSV 1, HSV 2, VZ and CMV, read using fluorescence microscope with the FITC/Evans Blue filter and the 40x objective.
- viii. For respiratory viruses, read using fluorescence microscope with the FITC/Evans Blue filter.

### Interpretation of Results

Positive: Bartel CMV monoclonal antibody: Bright apple green fluorescence of cytoplasmic inclusion (late antigen) and homogenous early nuclear antigen in CMV-CPE cells.

Chemicon SimulFluor Respiratory Screen:

All respiratory viruses except RSV show bright apple green fluorescence of the cytoplasm and/or nucleus of the infected cell.

RSV shows bright gold fluorescence of the cytoplasm and/or nucleus of the infected cell.

Chemicon SimulFluor Flu A/Flu B:

Influenzae A virus shows bright apple green fluorescence.

Influenzae B virus shows bright gold fluorescence.

Chemicon SimulFluor RSV/Para 3:

RSV virus shows bright apple green fluorescence.

Parainfluenzae 3 shows bright gold fluorescence.

Chemicon SimulFluor Para 123/Adeno:

Parainfluenza 1,2,3 viruses show bright apple green fluorescence.

Adenovirus shows bright gold fluorescence.

Chemicon individual monoclonal antibodies:

Parainfluenzae 1 and 2, and adenovirus show bright apple green fluorescence.

Negative: Red Cells with no apple-green fluorescence.

TML/MSH Microbiology Department Policy & Procedure Manual	<b>Policy # MI/VIR/16/05/v02</b>	Page 4 of 4
<b>Virology Manual</b>		

**IV. Quality Control**

Appropriate positive and negative control slides should be stained with each batch.

**VI. Reference**

Isenberg, H.D., 1992, ASM. Clinical Microbiology Procedures Handbook Vol. 2.