

TML/MSH Microbiology Department Policy & Procedure Manual	<b>Policy # MI/EDU/04/03/v01</b>	Page 1 of 1
Section: <b>Education Manual</b>	Subject Title: <b>Appendix III - Gram Stain Reactions</b>	
Issued by: <b>LABORATORY MANAGER</b>	Original Date: June 18, 2001	
Approved by: Laboratory Director	Revision Date:	

### **APPENDIX III - GRAM STAIN REACTIONS**

- |                        |   |
|------------------------|---|
| Gram positive cocci    | <ul style="list-style-type: none"> <li>- <i>Staphylococcus aureus</i></li> <li>- Coagulase negative <i>Staphylococcus</i></li> <li>- <i>Enterococcus faecalis</i></li> <li>- Viridans streptococcus</li> <li>- Beta-haemolytic group B streptococcus</li> </ul> |
| Gram positive bacilli  | <ul style="list-style-type: none"> <li>- <i>Clostridium novyi</i></li> <li>- <i>Bacillus</i> species</li> <li>- <i>Corynebacterium</i> species</li> </ul>   |
| Gram negative coccus   | <ul style="list-style-type: none"> <li>- <i>Moraxella catarhalis</i></li> <li>- <i>Neisseria</i> species</li> </ul>   |
| Gram negative bacillus | <ul style="list-style-type: none"> <li>- <i>Escherichia coli</i></li> <li>- <i>Proteus mirabilis</i></li> <li>- <i>Pseudomonas aeruginosa</i></li> </ul>  |